

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Research® plus Pipettes

Eppendorf Research® plus/ Eppendorf Research® plus Move It®	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader®/ Microloader™ 0.5-20 µL
epT.I.P.S.®	■		■	■	■				■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention®	■			■					■		■			■		■			
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■	■	■	■	■	■	
ep Dualfilter T.I.P.S.® LoRetention®	■			■				■		■	■								
ep Dualfilter T.I.P.S.® SealMax®	■			■				■	■	■	■								
Fixed volume, single-channel pipettes																			
10 µL	•	■	■	•															•
20 µL			•	■															•
10 µL							■	■	■	•									
20 µL	■						■	■	■	•									
25 µL/50 µL/100 µL	■							■	■	•									
200 µL	■								■	•									
200 µL/250 µL/500 µL/1,000 µL	■										■	•	•						
Variable volume, single-, 8- and 12-channel pipettes																			
0.1 – 2.5 µL	■	■	•	•															
0.5 – 10 µL	■	•	■	•															•
2 – 20 µL	■	◆	◆	■															•
2 – 20 µL	■						■	■	■	•									
10 – 100 µL	■						◆	■	■	•									
20 – 200 µL	■						◆	◆	■	•									
30 – 300 µL	■						◆	◆	◆	■									
100 – 1,000 µL	■										■	•	•						
120 – 1,200 µL	■											•	■						
0.25 – 2.5 mL	■													■					
0.5 – 5 mL	■														■	•			
1 – 10 mL	■																■	•	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
Move It®, variable volume, 4-, 6-, 8- and 12-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
30 – 300 µL	■						◆	◆	◆	■									
120 – 1,200 µL	■											•	■						

■ Test tip (other tips may require adjustment) • Full volume range ◆ Limited volume range

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)



For more information go to  
[www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Reference® 2 Pipettes

Eppendorf Reference® 2	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader®/ Microloader™ 0.5 – 20 µL
epT.I.P.S.®	■		■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention®	■			■			■		■								
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■		■	■	■	■	■		
ep Dualfilter T.I.P.S.® LoRetention®	■			■		■		■	■								
ep Dualfilter T.I.P.S.® SealMax®	■			■		■	■	■	■								
Fixed volume, single-channel pipettes																	
1 µL/2 µL	■	●	●														
5 µL/10 µL	■	■	■	●													●
20 µL	■		●	■													●
10 µL	■				■	■	■	●									
20 µL	■				■	■	■	●									
25 µL/50 µL/100 µL	■					■	■	●									
200 µL	■						■	●									
200 µL/250 µL/500 µL/1,000 µL	■								■	●	●						
2 mL/2.5 mL	■											■					
Variable volume, single-, 8- and 12-channel pipettes																	
0.1 – 2.5 µL	■	●	●														●
0.5 – 10 µL	■	■	■	●													●
2 – 20 µL	■	◆	●	■													●
2 – 20 µL	■				◆	◆	■	●									
10 – 100 µL	■				◆	■	■	●									
20 – 200 µL	■				■	■	■	●									
30 – 300 µL	■				◆	◆	◆	■									
100 – 1,000 µL	■								■	●	●						
0.25 – 2.5 mL	■											■					
0.5 – 5 mL	■												■	●			
1 – 10 mL	■														■	●	

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)



For more information go to  
[www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)

Combination of epT.I.P.S.® Pipette Tips with Eppendorf Xplorer®/Xplorer® plus Pipettes

Eppendorf Xplorer®/Eppendorf Xplorer® plus/ Eppendorf Xplorer® plus Move It®	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader®/ Microloader™ 0.5-20 µL
epT.I.P.S.®	■		■	■	■				■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S.® LoRetention®	■			■					■		■								
ep Dualfilter T.I.P.S.®	■	■		■	■	■	■	■	■	■	■		■	■	■	■	■	■	
ep Dualfilter T.I.P.S.® LoRetention®	■			■				■		■	■								
ep Dualfilter T.I.P.S.® SealMax®	■			■				■	■	■	■								
Variable volume, single-, 8- and 12-channel pipettes																			
0.5 – 10 µL	■	●	■	●															●
1 – 20 µL	■	◆	●	■															●
5 – 100 µL	■						◆	■	■	●									
10 – 200 µL	■						◆	◆	■	●									
15 – 300 µL	■						◆	◆	◆	■									
50 – 1,000 µL	■										■	●	●						
50 – 1,200 µL	■											■	●						
0.1 – 2.5 mL	■													■					
0.2 – 5 mL	■														■	●			
0.5 – 10 mL	■																■	●	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
Move It®, variable volume, 4-, 6-, 8- and 12-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
15 – 300 µL	■						◆	◆	◆	■									
50 – 1,200 µL	■											■	●						

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range

Your local distributor: [www.eppendorf.com/contact](http://www.eppendorf.com/contact)  
Eppendorf SE · 22339 Hamburg · Germany  
[eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com) · [www.eppendorf.com](http://www.eppendorf.com)

[www.eppendorf.com](http://www.eppendorf.com)



For more information go to  
[www.eppendorf.com/eptips](http://www.eppendorf.com/eptips)