

Lisa xxx

SOS serverist SSE sõnumi näidis sündmuse valimise kohta

```
{
  "SOS_EVENT": [
    {
      "tabEnum": "RESSURSS",
      "timestamp": 1723209891301,
      "status": "OPENED",
      "sosEventId": "139446",
      "sosEventVersion": 16,
      "sosEventNo": "2408080001",
      "aikEventId": null,
      "ambulance": {
        "active": true,
        "event": {
          "code": "32*",
          "name": "LIIKLUSÖNNETUS",
          "eventEnum":
"AIK_EVENT_AMBULANCE_TRAFFIC_ACCIDENT",
          "editorName": "Nipi Tiri",
          "editorIdCode": "47001010001",
          "lastUpdated": "2024-08-09T13:24:51+0000"
        },
        "priority": {
          "code": "PRIORITY_C",
          "recommendedCode": "PRIORITY_C",
          "editorName": "Nipi Tiri",
          "editorIdCode": "47001010001",
          "lastUpdated": "2024-08-09T13:24:51+0000"
        },
        "serviceStatus": "STATUS_OPEN_WITHOUT_RESOURCES",
        "criterias": [
          {
            "active": true,
            "criteriaId": "b513aa16-8eb0-4629-9a64-
2ca08304636d",
            "criteriaAnswerId": "a2251556-c7cf-
47e4-9736-fb4123eeb4d8",
            "originCriteriaId": null,
            "originAnswerId": null,
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-
08T07:05:26+0000"
          }
        ],
        "eventHist": [],
        "priorityHist": []
      },
      "rescue": {
        "active": true,
        "event": {
          "code": "16F",
          "name": "LIIKLUSÖNNETUS",
          "eventEnum":
"AIK_EVENT_RESCUE_TRAFFIC_ACCIDENT",
          "editorName": "Nipi Tiri",
          "editorIdCode": "47001010001",
```

```

        "lastUpdated": "2024-08-09T13:24:51+0000"
    },
    "priority": {
        "code": "PRIORITY_1",
        "recommendedCode": "PRIORITY_1",
        "editorName": "Nipi Tiri",
        "editorIdCode": "47001010001",
        "lastUpdated": "2024-08-09T13:24:51+0000"
    },
    "serviceStatus": "STATUS_OPEN_WITHOUT_RESOURCES",
    "criterias": [
        {
            "active": true,
            "criteriaId": "095c8823-0397-4626-934e-
929eefd7cd34",
            "criteriaAnswerId": null,
            "originCriteriaId": null,
            "originAnswerId": null,
            "editorName": null,
            "editorIdCode": null,
            "lastUpdated": null
        },
        {
            "active": false,
            "criteriaId": "f2491ccb-9038-4ald-9ea7-
bc4d11407477",
            "criteriaAnswerId": "c8165f70-ffe7-
4986-afea-5ca08e7b9beb",
            "originCriteriaId": null,
            "originAnswerId": null,
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-
09T13:23:59+0000"
        },
        {
            "active": true,
            "criteriaId": "507474a1-92f1-4f40-921f-
93ec83d36359",
            "criteriaAnswerId": "b8c23fa8-d346-
4e01-9845-b51211edc2ae",
            "originCriteriaId": null,
            "originAnswerId": null,
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-
09T13:23:57+0000"
        }
    ],
    "eventHist": [
        {
            "code": "20A",
            "name": "ABI OSUTAMINE",
            "eventEnum":
"AIK_EVENT_RESCUE_PROVISION_OF_AID",
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",

```

```

        "lastUpdated": "2024-08-
09T13:23:36+0000"
    },
    "priorityHist": [
        {
            "code": "PRIORITY_1A",
            "recommendedCode": "PRIORITY_1A",
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-
09T13:23:36+0000"
        }
    ],
    "police": {
        "active": true,
        "event": {
            "code": "B2",
            "name": "LIIKLUSÖNNETUS
VIGASTATUTE/HUKKUNUTEGA",
            "eventEnum":
"AIK_EVENT_POLICE_TRAFFIC_ACCIDENT_WITH_INJURED_DEAD_PEOPLE",
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-09T13:24:51+0000"
        },
        "priority": {
            "code": "PRIORITY_B",
            "recommendedCode": "PRIORITY_B",
            "editorName": "Nipi Tiri",
            "editorIdCode": "47001010001",
            "lastUpdated": "2024-08-09T13:24:51+0000"
        },
        "serviceStatus": "STATUS_OPEN_WITHOUT_RESOURCES",
        "criteria": [
            {
                "active": true,
                "criteriaId": "70835929-3d6e-4734-b08f-
fa6926dc392c",
                "criteriaAnswerId": null,
                "originCriteriaId": "b513aa16-8eb0-
4629-9a64-2ca08304636d",
                "originAnswerId": null,
                "editorName": "Nipi Tiri",
                "editorIdCode": "47001010001",
                "lastUpdated": "2024-08-
08T07:05:13+0000"
            }
        ],
        "eventHist": [],
        "priorityHist": []
    },
    "editorName": "Nipi Tiri",
    "editorIdCode": "47001010001",
    "lastUpdated": "2024-08-09T13:24:35+0000"
}
]

```

```
}
```

SOS serverist SSE sõnumi näidis sündmuse valiku eemaldamise kohta

```
{
  "SOS_EVENT": [
    {
      "tabEnum": "RESSURSS",
      "timestamp": 1723210311387,
      "status": "CLOSED",
      "sosEventId": 139446
    }
  ]
}

{
  "SOS_EVENT": [
    {
      "tabEnum": "RESSURSS",
      "timestamp": 1723210396695,
      "status": "CLOSED",
      "sosEventId": null
    }
  ]
}
```

Kiirabi ressursside registri REST teenuste spetsifikatsioon

openapi: 3.0.1

info:

title: Kiirabi ressursid API

version: 0.1.0

servers:

- url: http://localhost:8079

- url: https://kares.hk.dev

- url: https://kares.hk.test

- url: https://kares.hk.sise

tags:

- name: enums

description: Enumid

- name: hospital

description: Haiglatega seotud teenused

- name: userSetting

description: Kasutaja seadetega seotud teenused

- name: resource

description: Ressurssidega seotud teenused

- name: dispatchPlan

description: Väljasõiduplaaniga seotud teenused

paths:

/rest/v1/dispatchPlan:

get:

tags:

- dispatchPlan

summary: Väljasõiduplaani küsimine. Piiratud autenditud kasutajale (smitTokenId'ga).

description: Väljasõiduplaani küsimine. Piiratud autenditud kasutajale (smitTokenId'ga).

operationId: dispatchPlan/get/index

parameters:

```

- name: eventNo
  in: query
  description: Sündmuse number
  required: true
  schema:
    type: string
- name: Accept
  in: header
  description: Peab olema `application/json`
  required: true
  schema:
    type: string
    default: application/json
  example: application/json
- name: smitTokenID
  in: header
  description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId`'d
  required: true
  schema:
    type: string
  example: db0334e-538a-4f7f-b7a5-3e3d88da099b
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/CommonResponse'
        examples:
          dispatchPlan/get/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/get/index:
            value: |-
              {
                "dispatchPlanId" : 13,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : "19ec9563-0eaf-45ad-b11a-
c66ab4e2ed54",
                "dispatchRuleName" : "dispatchRuleName",
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ {
                  "active" : true,
                  "cancelled" : false,
                  "alarmed" : true,

```

```

        "done" : false,
        "extraResourceVehicleTypeEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "extraResourceVehicleTypeEnumText" :
"Arstibrigaad",
        "resource" : {
            "id" : 43,
            "externalId" : null,
            "name" : "NAME 2",
            "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
            "vehicleTypeAmbulanceEnumText" :
"Arstibrigaad",
            "stationId" : "STATION_ID",
            "distanceMeters" : null,
            "distanceSeconds" : null,
            "lastTimeOnEvent" : null
        }
    }, {
        "active" : true,
        "cancelled" : false,
        "alarmed" : true,
        "done" : false,
        "ruleResourceNeed" : {
            "orderNo" : 5,
            "amount" : 1,
            "conditionEnum" : "CONDITION_OR",
            "ruleResourceNeedTerms" : [ {
                "orderNo" : 1,
                "terms" : {
                    "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE"
                }
            }, {
                "orderNo" : 2,
                "terms" : {
                    "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
                }
            }
        ]
    },
        "extraResourceVehicleTypeEnum" : null,
        "extraResourceVehicleTypeEnumText" : null,
        "resource" : {
            "id" : 42,
            "externalId" : null,
            "name" : "NAME 1",
            "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
            "vehicleTypeAmbulanceEnumText" : "Öebrigaad",
            "stationId" : "STATION_ID",
            "distanceMeters" : 20.0,
            "distanceSeconds" : 20.0,
            "lastTimeOnEvent" : null
        }
    }
],
        "alternativeResources" : [ {
            "id" : 13,

```



```

        "DEVICE_TYPE_MOBILE_CHIEF" : "MOBIIL autovanem",
        "DEVICE_TYPE_NUGIS" : "NUGIS",
        "DEVICE_TYPE_TETRA" : "TETRA",
        "DEVICE_TYPE_TETRA_VEHICLE" : "TETRA autojaam",
        "DEVICE_TYPE_TETRA_CHIEF" : "TETRA autovanem",
        "DEVICE_TYPE_TETRA_PERSON" : "TETRA käsijaam",
        "DEVICE_TYPE_DERIVED" : "Tuletatud"
    },
    "resourceStatusEnums" : {
        "RESOURCE_STATUS_ANNULLED" : "Annulleeritud",
        "RESOURCE_STATUS_AWAY_BUSY" : "Eemal hõivatud",
        "RESOURCE_STATUS_AWAY" : "Eemal vaba",
        "RESOURCE_STATUS_AT_HOSPITAL" : "Haiglas",
        "RESOURCE_STATUS_TO_HOSPITAL" : "Haiglasse",
        "RESOURCE_STATUS_AT_HOME" : "Kodus",
        "RESOURCE_STATUS_ON_SITE" : "Kohal",
        "RESOURCE_STATUS_LEAVING_BUSY" : "Lahkun
hõivatud",

        "RESOURCE_STATUS_LEAVING" : "Lahkun vaba",
        "RESOURCE_STATUS_ALARMED" : "Teatamine",
        "RESOURCE_STATUS_ALARMED_BY_RADIO" : "Teatamine
raadio teel",

        "RESOURCE_STATUS_ACCEPTED" : "Vastuvõetud",
        "RESOURCE_STATUS_DEPARTURE" : "Väljasõit",
        "RESOURCE_STATUS_DEPARTURE_BY_RADIO" : "Väljasõit
raadio teel"
    },
    "eventStatusEnums" : {
        "STATUS_DELETED" : "Kustutatud",
        "STATUS_DRAFT" : "Mustand",
        "STATUS_OPEN" : "Töös",
        "STATUS_OPEN_WITHOUT_RESOURCES" : "Töös
ressursita",

        "STATUS_OPEN_PARTLY" : "Töös osaliselt",
        "STATUS_OPEN_WITH_RESOURCES" : "Töös ressursiga",
        "STATUS_OPEN_RESOURCES_LEFT" : "Töös ressursid
lahkunud",

        "STATUS_CLOSED" : "Lõpetatud",
        "STATUS_ANNULLED" : "Annulleeritud"
    },
    "priorityEnums" : {
        "PRIORITY_1" : {
            "value" : "1",
            "shortValue" : "1",
            "orderNo" : 100,
            "weight" : 100,
            "department" : [ "SERVICE_RESCUE" ]
        },
        "PRIORITY_1A" : {
            "value" : "1A",
            "shortValue" : "1A",
            "orderNo" : 90,
            "weight" : 90,
            "department" : [ "SERVICE_RESCUE" ]
        },
        "PRIORITY_2" : {
            "value" : "2",
            "shortValue" : "2",

```



```

        "orderNo" : 110,
        "weight" : 110,
        "department" : [ "SERVICE_RESCUE" ]
    },
    "PRIORITY_3" : {
        "value" : "3",
        "shortValue" : "3",
        "orderNo" : 120,
        "weight" : 120,
        "department" : [ "SERVICE_RESCUE" ]
    },
    "PRIORITY_4" : {
        "value" : "4",
        "shortValue" : "4",
        "orderNo" : 130,
        "weight" : 130,
        "department" : [ "SERVICE_RESCUE" ]
    },
    "PRIORITY_A" : {
        "value" : "A",
        "shortValue" : "A",
        "orderNo" : 140,
        "weight" : 140,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_B" : {
        "value" : "B",
        "shortValue" : "B",
        "orderNo" : 150,
        "weight" : 150,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_C" : {
        "value" : "C",
        "shortValue" : "C",
        "orderNo" : 160,
        "weight" : 160,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_D" : {
        "value" : "D",
        "shortValue" : "D",
        "orderNo" : 170,
        "weight" : 170,
        "department" : [ "SERVICE_AMBULANCE" ]
    },
    "PRIORITY_INFO" : {
        "value" : "Info",
        "shortValue" : "I",
        "orderNo" : 20,
        "weight" : 20,
        "department" : [ "SERVICE_POLICE" ]
    },
    "PRIORITY_CONSULTATION" : {
        "value" : "Nöuanne",

```

```

        "shortValue" : "N",
        "orderNo" : 10,
        "weight" : 10,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE", "SERVICE_RESCUE" ]
    }
},
"areaEnums" : {
    "AREA_EAST" : "Ida",
    "AREA_WEST" : "Lääne",
    "AREA_SOUTH" : "Lõuna",
    "AREA_NORTH" : "Põhja"
},
"conditionEnums" : {
    "CONDITION_AND" : "ja",
    "CONDITION_OR" : "või"
}
}
}
/rest/v1/hospital:
get:
    tags:
    - hospital
    summary: Haiglata pearing. Piiratud autenditud kasutajale
(smitTokenId'ga).
    description: Haiglata pearing. Piiratud autenditud kasutajale
(smitTokenId'ga).
    operationId: hospital/list
    parameters:
    - name: Accept
      in: header
      description: Peab olema `application/json`
      required: true
      schema:
        type: string
        default: application/json
      example: application/json
    - name: smitTokenID
      in: header
      description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
      `smitTokenId`'d
      required: true
      schema:
        type: string
      example: aldf644-0b99-4aae-99f1-fadb75ff3f80
    responses:
      "200":
        description: "200"
        content:
          application/json;charset=UTF-8:
            schema:
              $ref: '#/components/schemas/HospitalResponse'
            examples:
              hospital/list:
                value: |-
                  {
                    "total" : 10,
                    "list" : [ {

```

```
"id" : 1,
"version" : 0,
"name" : "HOSPITAL_0",
"code" : "CODE_0",
"address" : null,
"addressWgsLongitude" : null,
"addressWgsLatitude" : null,
"areaEnum" : null,
"areaEnumText" : null,
"description" : null,
"deleted" : false,
"creatorSessionName" : "Nipi Tiri",
"creatorSessionIdCode" : "47001010001",
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"dateCreated" : "2024-08-08T17:02:32+03",
"lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 2,
  "version" : 0,
  "name" : "HOSPITAL_1",
  "code" : "CODE_1",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 3,
  "version" : 0,
  "name" : "HOSPITAL_2",
  "code" : "CODE_2",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 4,
  "version" : 0,
  "name" : "HOSPITAL_3",
  "code" : "CODE_3",
```

```
"address" : null,
"addressWgsLongitude" : null,
"addressWgsLatitude" : null,
"areaEnum" : null,
"areaEnumText" : null,
"description" : null,
"deleted" : false,
"creatorSessionName" : "Nipi Tiri",
"creatorSessionIdCode" : "47001010001",
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"dateCreated" : "2024-08-08T17:02:32+03",
"lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 5,
  "version" : 0,
  "name" : "HOSPITAL_4",
  "code" : "CODE_4",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 6,
  "version" : 0,
  "name" : "HOSPITAL_5",
  "code" : "CODE_5",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 7,
  "version" : 0,
  "name" : "HOSPITAL_6",
  "code" : "CODE_6",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
```

```
"areaEnumText" : null,
"description" : null,
"deleted" : false,
"creatorSessionName" : "Nipi Tiri",
"creatorSessionIdCode" : "47001010001",
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"dateCreated" : "2024-08-08T17:02:32+03",
"lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 8,
  "version" : 0,
  "name" : "HOSPITAL_7",
  "code" : "CODE_7",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 9,
  "version" : 0,
  "name" : "HOSPITAL_8",
  "code" : "CODE_8",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:32+03",
  "lastUpdated" : "2024-08-08T17:02:32+03"
}, {
  "id" : 10,
  "version" : 0,
  "name" : "HOSPITAL_9",
  "code" : "CODE_9",
  "address" : null,
  "addressWgsLongitude" : null,
  "addressWgsLatitude" : null,
  "areaEnum" : null,
  "areaEnumText" : null,
  "description" : null,
  "deleted" : false,
  "creatorSessionName" : "Nipi Tiri",
```

```

        "creatorSessionIdCode" : "47001010001",
        "editorSessionName" : "Nipi Tiri",
        "editorSessionIdCode" : "47001010001",
        "dateCreated" : "2024-08-08T17:02:32+03",
        "lastUpdated" : "2024-08-08T17:02:32+03"
    } ]
}
/rest/v1/resource:
get:
  tags:
  - resource
  summary: Ressurside päring (aktiivsete). Piiratud autenditud
kasutajale (smitTokenId'ga).
  description: Ressurside päring (aktiivsete). Piiratud autenditud
kasutajale
(smitTokenId'ga).
  operationId: resource/list
  parameters:
  - name: id
    in: query
    description: Ressursi ID
    required: false
    schema:
      type: string
  - name: limit
    in: query
    description: "Lehekülgedeks jagamiseks, ühe lehekülje kirjete
arv. Number,\
    \ vaike väärtus on 10 000"
    required: false
    schema:
      type: string
  - name: offset
    in: query
    description: "Lehekülgedeks jagamiseks, mitmendast kirjest
alates. Number,\
    \ vaikumisi 0"
    required: false
    schema:
      type: string
  - name: sort
    in: query
    description: "Komadega eraldatud nimekiri väljade nimedest,
milledel alusel\
    \ tuleb tulemus sorteerida. Vaikumisi `id`\n\rNt.
`sort=name,firstName`\
    \ sorteerib tulemuse väljade `name` ja seejärel `firstName`
järgi."
    required: false
    schema:
      type: string
  - name: order
    in: query
    description: "Komadega eraldatud nimekiri, mis suunas eelnevalt
antud välju\
    \ sorteerida - asc/desc (kasvavalt/kahanevalt). Vaikumisi
`asc`\n\rNt. `sort=name,firstName&order=asc,desc`\n\

```

```

    \rTähendab, et esmalt sorteeritakse välja `name` kasvavalt ja
seejärel vä\
    lja `firstName` kahanevas järjekorras."
required: false
schema:
  type: string
  enum:
    - asc
    - desc
- name: filter
  in: query
  description: "Või (OR) filtri tingimused. Kõiki välju, mis on
teenuse juures\
  \ filtri tingimustena kirjeldatud on võimalik kasutada ka
omavahel kui VÕ\
  I (OR) tingimusena.\n\rTingimused on teineteisest eraldatud
komaga (mitte\
  \ `&` märgiga) ja väärtus järgneb koolonile (:) mitte võrdus
märgile (=).\n\
  \rNt.
`filter=showToEveryone:true,showToUser:true,notificationTime:gt:2017-11-
01T12:00:00+03,notificationTime:lt:2017-11-01T12:00:00+03`\n\
  \rTähendab - (\n\r 'showToEveryone' = true VÕI \n\r
'showToUser' = true\
  \ VÕI \n\r 'notificationTime' > '2017-11-01T12:00:00+03' VÕI
\n\r 'notificationTime'\
  \ < '2017-11-01T12:00:00+03'\n\r)"
required: false
schema:
  type: string
- name: Accept
  in: header
  description: Peab olema `application/json`
required: true
schema:
  type: string
  default: application/json
  example: application/json
- name: smitTokenID
  in: header
  description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId`'d
required: true
schema:
  type: string
  example: 89c93eec-8a24-4fb7-aa58-7ffe69d64776
responses:
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/ResourceResponse'
        examples:
          resource/list:
            value: |-
              {

```

```
"total" : 2,
"list" : [ {
  "id" : 19,
  "version" : 3,
  "name" : "RESOURCE_NAME_1",
  "active" : true,
  "station" : "station",
  "stationId" : "stationId",
  "stationLongitude" : null,
  "stationLatitude" : null,
  "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE",
  "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
  "areaEnum" : "AREA_WEST",
  "areaEnumText" : "Lääne",
  "responseRange" : null,
  "responseTime" : null,
  "externalId" : null,
  "comment" : null,
  "expiring" : false,
  "validUntil" : null,
  "sosVersion" : null,
  "resourceLocation" : {
    "id" : 19,
    "version" : 0,
    "deviceTypeEnum" : "DEVICE_TYPE_TETRA",
    "deviceTypeEnumText" : "TETRA",
    "locationTime" : null,
    "longitude" : null,
    "latitude" : null,
    "speed" : null,
    "direction" : null,
    "lastUpdated" : "2024-08-08T17:02:35+03"
  },
  "resourceStatus" : {
    "id" : 19,
    "version" : 0,
    "resourceStatusEnum" :
"RESOURCE_STATUS_DEPARTURE",
    "resourceStatusEnumText" : "Väljasõit",
    "statusNote" : null,
    "lastEventDate" : null,
    "editorSessionName" : "Nipi Tiri",
    "editorSessionIdCode" : "47001010001",
    "lastUpdated" : "2024-08-08T17:02:35+03"
  },
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03",
  "devices" : [ {
    "id" : 63,
    "version" : 0,
    "resourceId" : 19,
    "resourceMemberId" : null,
    "externalId" : "RESOURCE_NAME_1",
```



```
"deviceTypeEnum" : "DEVICE_TYPE_MGIS",
"deviceTypeEnumText" : "MGIS",
"deviceLocation" : {
  "id" : 9,
  "version" : 0,
  "locationTime" : "2022-10-17T16:56:00+03",
  "longitude" : 25.363565,
  "latitude" : 58.456746,
  "speed" : 5.0,
  "direction" : 33.0,
  "lastUpdated" : "2024-08-08T17:02:35+03"
},
"dateCreated" : "2024-08-08T17:02:35+03",
"lastUpdated" : "2024-08-08T17:02:35+03"
}, {
  "id" : 61,
  "version" : 0,
  "resourceId" : 19,
  "resourceMemberId" : 1,
  "externalId" : "2500001",
  "deviceTypeEnum" : "DEVICE_TYPE_TETRA",
  "deviceTypeEnumText" : "TETRA",
  "deviceLocation" : {
    "id" : 7,
    "version" : 0,
    "locationTime" : null,
    "longitude" : null,
    "latitude" : null,
    "speed" : null,
    "direction" : null,
    "lastUpdated" : "2024-08-08T17:02:35+03"
  },
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
}, {
  "id" : 62,
  "version" : 0,
  "resourceId" : 19,
  "resourceMemberId" : 1,
  "externalId" : "email@email.fake",
  "deviceTypeEnum" : "DEVICE_TYPE_EMAIL",
  "deviceTypeEnumText" : "E-mail",
  "deviceLocation" : {
    "id" : 8,
    "version" : 0,
    "locationTime" : null,
    "longitude" : null,
    "latitude" : null,
    "speed" : null,
    "direction" : null,
    "lastUpdated" : "2024-08-08T17:02:35+03"
  },
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
} ],
"events" : [ {
  "id" : 4,
  "version" : 0,
```

```
"eventId" : 9,
"resourceId" : 19,
"radio" : false,
"alarmedDate" : null,
"acceptedDate" : null,
"departureDate" : "2024-08-08T17:02:35+03",
"onSiteDate" : null,
"hospitalDate" : null,
"hospitalId" : null,
"hospitalName" : null,
"leavingDate" : null,
"homeDate" : null,
"annulledDate" : null,
"additionalData" : null,
"stationAlarmedStartDate" : null,
"stationAlarmedEndDate" : null,
"estimatedOnSiteArrivalTimeInSeconds" : null,
"active" : true,
"externalId" : null,
"removedFromSync" : false,
"statusEnum" : "RESOURCE_STATUS_DEPARTURE",
"statusEnumText" : "Väljasõit",
"creatorSessionName" : "Nipi Tiri",
"creatorSessionIdCode" : "47001010001",
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"dateCreated" : "2024-08-08T17:02:35+03",
"lastUpdated" : "2024-08-08T17:02:35+03"
} ],
"members" : [ {
  "id" : 1,
  "version" : 2,
  "resourceId" : 19,
  "externalId" : "1158545",
  "isChief" : true,
  "active" : true,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
}, {
  "id" : 2,
  "version" : 0,
  "resourceId" : 19,
  "externalId" : "1158546",
  "isChief" : false,
  "active" : true,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
} ]
}, {
  "id" : 20,
```

```

"version" : 1,
"name" : "RESOURCE_NAME_2",
"active" : true,
"station" : "station",
"stationId" : "stationId",
"stationLongitude" : null,
"stationLatitude" : null,
"vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE",
"vehicleTypeAmbulanceEnumText" : "Õebrigaad",
"areaEnum" : "AREA_WEST",
"areaEnumText" : "Lääne",
"responseRange" : null,
"responseTime" : null,
"externalId" : null,
"comment" : null,
"expiring" : false,
"validUntil" : null,
"sosVersion" : null,
"resourceLocation" : {
  "id" : 20,
  "version" : 0,
  "deviceTypeEnum" : "DEVICE_TYPE_TETRA",
  "deviceTypeEnumText" : "TETRA",
  "locationTime" : null,
  "longitude" : null,
  "latitude" : null,
  "speed" : null,
  "direction" : null,
  "lastUpdated" : "2024-08-08T17:02:35+03"
},
"resourceStatus" : {
  "id" : 20,
  "version" : 0,
  "resourceStatusEnum" :
"RESOURCE_STATUS_TO_HOSPITAL",
"resourceStatusEnumText" : "Haiglasse",
"statusNote" : null,
"lastEventDate" : null,
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"lastUpdated" : "2024-08-08T17:02:35+03"
},
"creatorSessionName" : "Nipi Tiri",
"creatorSessionIdCode" : "47001010001",
"editorSessionName" : "Nipi Tiri",
"editorSessionIdCode" : "47001010001",
"dateCreated" : "2024-08-08T17:02:35+03",
"lastUpdated" : "2024-08-08T17:02:35+03",
"devices" : [ {
  "id" : 64,
  "version" : 0,
  "resourceId" : 20,
  "resourceMemberId" : null,
  "externalId" : "RESOURCE_NAME_2",
  "deviceTypeEnum" : "DEVICE_TYPE_MGIS",
  "deviceTypeEnumText" : "MGIS",
  "deviceLocation" : {

```

```
    "id" : 10,
    "version" : 0,
    "locationTime" : "2022-10-17T16:56:00+03",
    "longitude" : 25.363565,
    "latitude" : 58.456746,
    "speed" : 5.0,
    "direction" : 33.0,
    "lastUpdated" : "2024-08-08T17:02:35+03"
  },
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
} ],
"events" : [ {
  "id" : 5,
  "version" : 0,
  "eventId" : 10,
  "resourceId" : 20,
  "radio" : false,
  "alarmedDate" : null,
  "acceptedDate" : null,
  "departureDate" : null,
  "onSiteDate" : null,
  "hospitalDate" : "2024-08-08T17:02:35+03",
  "hospitalId" : 21,
  "hospitalName" : "name",
  "leavingDate" : null,
  "homeDate" : null,
  "annulledDate" : null,
  "additionalData" : null,
  "stationAlarmedStartDate" : null,
  "stationAlarmedEndDate" : null,
  "estimatedOnSiteArrivalTimeInSeconds" : null,
  "active" : true,
  "externalId" : null,
  "removedFromSync" : false,
  "statusEnum" : "RESOURCE_STATUS_TO_HOSPITAL",
  "statusEnumText" : "Haiglassse",
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
} ],
"members" : [ {
  "id" : 3,
  "version" : 0,
  "resourceId" : 20,
  "externalId" : "1258545",
  "isChief" : true,
  "active" : true,
  "creatorSessionName" : "Nipi Tiri",
  "creatorSessionIdCode" : "47001010001",
  "editorSessionName" : "Nipi Tiri",
  "editorSessionIdCode" : "47001010001",
  "dateCreated" : "2024-08-08T17:02:35+03",
  "lastUpdated" : "2024-08-08T17:02:35+03"
} ]
```

```

        } ]
    }
    "401":
      description: "401"
      content:
        application/json;charset=UTF-8:
          schema:
            $ref: '#/components/schemas/CommonResponse'
          examples:
            resource/list-forbidden:
              value: |-
                {
                  "status" : 401,
                  "message" : "Kasutajal puudub kehtiv sessioon"
                }
/rest/v1/userSetting:
  get:
    tags:
      - userSetting
    summary: "Kasutajaseadete pärimine. Vajab autenditud kasutajat
(smitTokenId),\
  \ kelle seadeid tagastatakse"
    description: "Kasutajaseadete pärimine. Vajab autenditud kasutajat
(smitTokenId),\
  \ kelle seadeid tagastatakse"
    operationId: userSetting/index
    parameters:
      - name: typeEnum
        in: query
        description: Seadete tunnus
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
          example: application/json
      - name: smitTokenID
        in: header
        description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId`'d
        required: true
        schema:
          type: string
          example: e19679d4-24a2-4f4a-881f-d9bf8dfe3f63
    responses:
      "204":
        description: "204"
      "401":
        description: "401"
        content:
          application/json;charset=UTF-8:
            schema:

```

```

    $ref: '#/components/schemas/CommonResponse'
  examples:
    userSetting/index-forbidden:
      value: |-
        {
          "status" : 401,
          "message" : "Kasutajal puudub kehtiv sessioon"
        }
"200":
  description: "200"
  content:
    application/json;charset=utf-8:
      schema:
        $ref: '#/components/schemas/rest-v1-userSetting486549215'
      examples:
        userSetting/index:
          value: |-
            {
              "filter1" : "93d52191-a95a-4b96-acf8-357f64c52376",
              "filter2" : [ "9d2b8929-c318-4ab3-9752-
fb4e0b1f4216", "51c11812-1c39-4d7f-b609-0ebed5ca7676", "570b2f2d-3cc5-
4773-965e-e295e025eeb6" ]
            }
"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/CommonResponse'
      examples:
        userSetting/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "typeEnum pole lubatud väärtusega,
null"
            }
/rest/v1/dispatchPlan/acceptResourceNeed:
  post:
    tags:
      - dispatchPlan
    summary: Ressursi vajaduse aktsepteerimine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    description: Ressursi vajaduse aktsepteerimine. Piiratud autenditud
kasutajale
      (smitTokenId'ga).
    operationId: dispatchPlan/acceptResourceNeed/index
    parameters:
      - name: eventNo
        in: query
        description: Sündmuse number
        required: true
        schema:
          type: string
      - name: resourceNeedId
        in: query
        description: Ressursi vajaduse tunnus
        required: true

```

```

    schema:
      type: string
- name: Accept
  in: header
  description: Peab olema `application/json`
  required: true
  schema:
    type: string
    default: application/json
    example: application/json
- name: smitTokenID
  in: header
  description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId`'d
  required: true
  schema:
    type: string
    example: 938de407-1dbc-4d67-aae9-ca4e9c2d8577
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/acceptResourceNeed/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/acceptResourceNeed/index:
            value: |-
              {
                "dispatchPlanId" : 6,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : null,
                "dispatchRuleName" : null,
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ ],
                "alternativeResources" : [ ]
              }
  "412":
    description: "412"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'

```

```

examples:
  dispatchPlan/acceptResourceNeed/index-notValid:
    value: |-
      {
        "status" : 412,
        "message" : "\"resourceNeedId\" peab olema
määratud",
        "errorHeader" : "Ressursi vajaduse aktsepteerimine
ebaõnnestus!"
      }
/rest/v1/dispatchPlan/addResource:
  post:
    tags:
      - dispatchPlan
    summary: Ressursi lisamine vajadusele. Piiratud autenditud
kasutajale (smitTokenId'ga).
    description: Ressursi lisamine vajadusele. Piiratud autenditud
kasutajale (smitTokenId'ga).
    operationId: dispatchPlan/addResource/index
    parameters:
      - name: eventNo
        in: query
        description: Sündmuse number
        required: true
        schema:
          type: string
      - name: resourceId
        in: query
        description: Ressursi tunnus
        required: false
        schema:
          type: string
      - name: resourceExternalId
        in: query
        description: Ressursi väline tunnus
        required: false
        schema:
          type: string
      - name: resourceName
        in: query
        description: Ressursi kutse
        required: false
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
          example: application/json
      - name: smitTokenID
        in: header
        description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
`smitTokenId`'d
        required: true

```



```

schema:
  type: string
  example: ed23b9a0-087f-4567-821e-9f458785b75e
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/addResource/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/addResource/index:
            value: |-
              {
                "dispatchPlanId" : 7,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : null,
                "dispatchRuleName" : null,
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ {
                  "active" : true,
                  "cancelled" : false,
                  "alarmed" : false,
                  "done" : false,
                  "extraResourceVehicleTypeEnum" : null,
                  "extraResourceVehicleTypeEnumText" : null,
                  "resource" : {
                    "id" : 30,
                    "externalId" : null,
                    "name" : "NAME",
                    "vehicleTypeEnum" :
"VEHICLE_TYPE_AMBULANCE",
                    "vehicleTypeEnumText" : "Reanimobiil",
                    "stationId" : "STATION_ID",
                    "distanceMeters" : 5408.8,
                    "distanceSeconds" : 292.8,
                    "lastTimeOnEvent" : null
                  }
                } ],
                "alternativeResources" : [ {
                  "id" : 5,
                  "externalId" : null,
                  "name" : "NAME",
                  "stationId" : "STATION_ID",

```

```

        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
        "distanceMeters" : 5408.8,
        "distanceSeconds" : 292.8,
        "lastTimeOnEvent" : null
    } ]
}
"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/addResource/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "Ressursi lisamiseks peab olema
määratud kas resourceId, resourceExternalId või resourceName!",
              "errorHeader" : "Ressursi lisamine ebaõnnestus!"
            }
/rest/v1/dispatchPlan/cancelResourceNeed:
  post:
    tags:
      - dispatchPlan
    summary: Ressursi vajaduse eemaldamine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    description: Ressursi vajaduse eemaldamine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    operationId: dispatchPlan/cancelResourceNeed/index
    parameters:
      - name: eventNo
        in: query
        description: Sündmuse number
        required: true
        schema:
          type: string
      - name: resourceNeedId
        in: query
        description: Ressursi vajaduse tunnus
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
          example: application/json
      - name: smitTokenID
        in: header
        description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
`smitTokenId`'d

```

```

required: true
schema:
  type: string
example: eb5a6ab8-15e7-45f1-b143-d5c20ca2b5e9
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/cancelResourceNeed/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/cancelResourceNeed/index:
            value: |-
              {
                "dispatchPlanId" : 10,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : null,
                "dispatchRuleName" : null,
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ {
                  "active" : true,
                  "cancelled" : false,
                  "alarmed" : false,
                  "done" : false,
                  "ruleResourceNeed" : {
                    "orderNo" : 2,
                    "amount" : 1,
                    "conditionEnum" : "CONDITION_AND",
                    "ruleResourceNeedTerms" : [ {
                      "orderNo" : 1,
                      "terms" : {
                        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
                      }
                    } ]
                  },
                  "extraResourceVehicleTypeEnum" : null,
                  "extraResourceVehicleTypeEnumText" : null,
                  "resource" : {
                    "id" : 34,
                    "externalId" : null,
                    "name" : "NAME 2",

```

```

        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
        "stationId" : "STATION_ID",
        "distanceMeters" : 5408.8,
        "distanceSeconds" : 592.8,
        "lastTimeOnEvent" : null
    }
} ],
"alternativeResources" : [ {
    "id" : 8,
    "externalId" : null,
    "name" : "NAME 3",
    "stationId" : "STATION_ID",
    "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Õebrigaad",
        "distanceMeters" : 25408.8,
        "distanceSeconds" : 1292.8,
        "lastTimeOnEvent" : null
    } ]
}
"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/cancelResourceNeed/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "\"resourceNeedId\" peab olema
määratud",
              "errorHeader" : "Ressursi vajaduse eemaldamine
ebaõnnestus!"
            }
/rest/v1/dispatchPlan/recalculate:
  post:
    tags:
      - dispatchPlan
    summary: Väljasõiduplaani ümberarvutamine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    description: Väljasõiduplaani ümberarvutamine. Piiratud autenditud
kasutajale
      (smitTokenId'ga).
    operationId: dispatchPlan/recalculate/index
    parameters:
      - name: eventNo
        in: query
        description: Sündmuse number
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`

```

```

    required: true
    schema:
      type: string
      default: application/json
      example: application/json
- name: smitTokenID
  in: header
  description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId`'d
  required: true
  schema:
    type: string
    example: 564deeb5-c323-4d43-bc15-483e74d03fc7
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/recalculate/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/recalculate/index:
            value: |-
              {
                "dispatchPlanId" : 16,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : null,
                "dispatchRuleName" : null,
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ {
                  "active" : true,
                  "cancelled" : false,
                  "alarmed" : false,
                  "done" : false,
                  "ruleResourceNeed" : {
                    "orderNo" : 1,
                    "amount" : 1,
                    "conditionEnum" : "CONDITION_AND",
                    "ruleResourceNeedTerms" : [ {
                      "orderNo" : 1,
                      "terms" : {
                        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"

```

```

    }
  } ]
},
"extraResourceVehicleTypeEnum" : null,
"extraResourceVehicleTypeEnumText" : null
}, {
  "active" : true,
  "cancelled" : false,
  "alarmed" : false,
  "done" : false,
  "ruleResourceNeed" : {
    "orderNo" : 2,
    "amount" : 1,
    "conditionEnum" : "CONDITION_AND",
    "ruleResourceNeedTerms" : [ {
      "orderNo" : 1,
      "terms" : {
        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_MEDICAL_BRIGADE"
      }
    } ]
  },
  "extraResourceVehicleTypeEnum" : null,
  "extraResourceVehicleTypeEnumText" : null
}, {
  "active" : true,
  "cancelled" : false,
  "alarmed" : false,
  "done" : false,
  "ruleResourceNeed" : {
    "orderNo" : 3,
    "amount" : 1,
    "conditionEnum" : "CONDITION_AND",
    "ruleResourceNeedTerms" : [ {
      "orderNo" : 1,
      "terms" : {
        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_MEDICAL_BRIGADE"
      }
    } ]
  },
  "extraResourceVehicleTypeEnum" : null,
  "extraResourceVehicleTypeEnumText" : null
}, {
  "active" : true,
  "cancelled" : false,
  "alarmed" : false,
  "done" : false,
  "ruleResourceNeed" : {
    "orderNo" : 4,
    "amount" : 1,
    "conditionEnum" : "CONDITION_OR",
    "ruleResourceNeedTerms" : [ {
      "orderNo" : 1,
      "terms" : {
        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
      }
    } ]
  },
  "extraResourceVehicleTypeEnum" : null,
  "extraResourceVehicleTypeEnumText" : null
}

```

```

        }, {
            "orderNo" : 2,
            "terms" : {
                "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE"
            }
        } ]
    },
    "extraResourceVehicleTypeEnum" : null,
    "extraResourceVehicleTypeEnumText" : null
}, {
    "active" : true,
    "cancelled" : false,
    "alarmed" : false,
    "done" : false,
    "ruleResourceNeed" : {
        "orderNo" : 5,
        "amount" : 1,
        "conditionEnum" : "CONDITION_OR",
        "ruleResourceNeedTerms" : [ {
            "orderNo" : 1,
            "terms" : {
                "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE"
            }
        }, {
            "orderNo" : 2,
            "terms" : {
                "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
            }
        } ]
    },
    "extraResourceVehicleTypeEnum" : null,
    "extraResourceVehicleTypeEnumText" : null
} ],
"alternativeResources" : [ {
    "id" : 16,
    "externalId" : null,
    "name" : "NAME",
    "stationId" : "STATION_ID",
    "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
    "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
    "distanceMeters" : 5408.8,
    "distanceSeconds" : 292.8,
    "lastTimeOnEvent" : null
} ]
}

```

```

"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/recalculate/index-notValid:
          value: |-

```

```

        {
            "status" : 412,
            "message" : "\"eventNo\" peab olema määratud",
            "errorHeader" : "Väljasõiduplaani värskendamine
ebaõnnestus!"
        }
    /rest/v1/dispatchPlan/removeResource:
    post:
        tags:
        - dispatchPlan
        summary: Ressursi eemaldamine vajaduselt. Piiratud autenditud
kasutajale (smitTokenId'ga).
        description: Ressursi eemaldamine vajaduselt. Piiratud autenditud
kasutajale
        (smitTokenId'ga).
        operationId: dispatchPlan/removeResource/index
        parameters:
        - name: eventNo
          in: query
          description: Sündmuse number
          required: true
          schema:
            type: string
        - name: resourceNeedId
          in: query
          description: Ressursi vajaduse tunnus
          required: true
          schema:
            type: string
        - name: Accept
          in: header
          description: Peab olema `application/json`
          required: true
          schema:
            type: string
            default: application/json
            example: application/json
        - name: smitTokenID
          in: header
          description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
          `smitTokenId`'d
          required: true
          schema:
            type: string
            example: 9cd4abce-9715-4332-9d32-14df1e2dd820
        responses:
            "401":
                description: "401"
                content:
                    application/json;charset=UTF-8:
                    schema:
                        $ref: '#/components/schemas/DispatchPlanResponse'
                    examples:
                        dispatchPlan/removeResource/index-forbidden:
                            value: |-
                                {
                                    "status" : 401,

```



```

        "message" : "Kasutajal puudub kehtiv sessioon"
    }
"200":
  description: "200"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/removeResource/index:
          value: |-
            {
              "dispatchPlanId" : 19,
              "eventNo" : "1",
              "eventVersion" : 0,
              "dispatchRuleId" : null,
              "dispatchRuleName" : null,
              "preferredResourceStationIds" : null,
              "resourceNeeds" : [ {
                "active" : true,
                "cancelled" : false,
                "alarmed" : false,
                "done" : false,
                "ruleResourceNeed" : {
                  "orderNo" : 1,
                  "amount" : 1,
                  "conditionEnum" : "CONDITION_OR",
                  "ruleResourceNeedTerms" : [ {
                    "orderNo" : 1,
                    "terms" : {
                      "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
                    }
                  }, {
                    "orderNo" : 2,
                    "terms" : {
                      "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE"
                    }
                  }
                ]
              },
              "extraResourceVehicleTypeEnum" : null,
              "extraResourceVehicleTypeEnumText" : null
            } ],
          "alternativeResources" : [ {
            "id" : 19,
            "externalId" : null,
            "name" : "NAME 2",
            "stationId" : "STATION_ID",
            "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
            "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
            "distanceMeters" : 25408.8,
            "distanceSeconds" : 1292.8,
            "lastTimeOnEvent" : null
          }, {
            "id" : 20,
            "externalId" : null,

```

```

        "name" : "NAME 3",
        "stationId" : "STATION_ID",
        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Õebrigaad",
        "distanceMeters" : 15408.8,
        "distanceSeconds" : 692.8,
        "lastTimeOnEvent" : null
    } ]
}
"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/removeResource/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "\"resourceNeedId\" peab olema
määratud",
              "errorHeader" : "Ressursi eemaldamine ebaõnnestus!"
            }
/rest/v1/dispatchPlan/restoreResourceNeed:
  post:
    tags:
      - dispatchPlan
    summary: Ressursi vajaduse taastamine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    description: Ressursi vajaduse taastamine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    operationId: dispatchPlan/restoreResourceNeed/index
    parameters:
      - name: eventNo
        in: query
        description: Sündmuse number
        required: true
        schema:
          type: string
      - name: resourceNeedId
        in: query
        description: Ressursi vajaduse tunnus
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
        example: application/json
      - name: smitTokenID
        in: header

```

description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali

```
  `smitTokenId` 'd
  required: true
  schema:
    type: string
  example: 56e40f84-3400-4eb0-b197-70f6744628ce
responses:
  "401":
    description: "401"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/restoreResourceNeed/index-forbidden:
            value: |-
              {
                "status" : 401,
                "message" : "Kasutajal puudub kehtiv sessioon"
              }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/DispatchPlanResponse'
        examples:
          dispatchPlan/restoreResourceNeed/index:
            value: |-
              {
                "dispatchPlanId" : 22,
                "eventNo" : "1",
                "eventVersion" : 0,
                "dispatchRuleId" : null,
                "dispatchRuleName" : null,
                "preferredResourceStationIds" : null,
                "resourceNeeds" : [ {
                  "active" : true,
                  "cancelled" : false,
                  "alarmed" : false,
                  "done" : false,
                  "ruleResourceNeed" : {
                    "orderNo" : 1,
                    "amount" : 1,
                    "conditionEnum" : "CONDITION_AND",
                    "ruleResourceNeedTerms" : [ {
                      "orderNo" : 1,
                      "terms" : {
                        "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
                      }
                    } ]
                  },
                  "extraResourceVehicleTypeEnum" : null,
                  "extraResourceVehicleTypeEnumText" : null,
                  "resource" : {
                    "id" : 65,
```

```

        "externalId" : null,
        "name" : "NAME 3",
        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
        "stationId" : "STATION_ID",
        "distanceMeters" : 15408.8,
        "distanceSeconds" : 692.8,
        "lastTimeOnEvent" : null
    }
}, {
    "active" : true,
    "cancelled" : false,
    "alarmed" : false,
    "done" : false,
    "ruleResourceNeed" : {
        "orderNo" : 2,
        "amount" : 1,
        "conditionEnum" : "CONDITION_OR",
        "ruleResourceNeedTerms" : [ {
            "orderNo" : 1,
            "terms" : {
                "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_REANIMOBILE"
            }
        } ], {
            "orderNo" : 2,
            "terms" : {
                "vehicleTypeCode" :
"VEHICLE_TYPE_AMBULANCE_SISTER_BRIGADE"
            }
        } ]
    },
    "extraResourceVehicleTypeEnum" : null,
    "extraResourceVehicleTypeEnumText" : null,
    "resource" : {
        "id" : 64,
        "externalId" : null,
        "name" : "NAME 2",
        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
        "stationId" : "STATION_ID",
        "distanceMeters" : 25408.8,
        "distanceSeconds" : 1292.8,
        "lastTimeOnEvent" : null
    }
} ],
"alternativeResources" : [ {
    "id" : 28,
    "externalId" : null,
    "name" : "NAME 2",
    "stationId" : "STATION_ID",
    "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
    "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
    "distanceMeters" : 25408.8,
    "distanceSeconds" : 1292.8,

```

```

        "lastTimeOnEvent" : null
    }, {
        "id" : 29,
        "externalId" : null,
        "name" : "NAME 3",
        "stationId" : "STATION_ID",
        "vehicleTypeAmbulanceEnum" :
"VEHICLE_TYPE_AMBULANCE",
        "vehicleTypeAmbulanceEnumText" : "Reanimobiil",
        "distanceMeters" : 15408.8,
        "distanceSeconds" : 692.8,
        "lastTimeOnEvent" : null
    } ]
}
"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/DispatchPlanResponse'
      examples:
        dispatchPlan/restoreResourceNeed/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "\"resourceNeedId\" peab olema
määratud",
              "errorHeader" : "Ressursi vajaduse taastamine
ebaõnnestus!"
            }
/rest/v1/enums/areaEnums:
  get:
    tags:
      - enums
    summary: Kõikide AreaEnum enumite päring
    description: Kõikide AreaEnum enumite päring
    operationId: enum/areaEnum
    parameters:
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
        example: application/json
    responses:
      "200":
        description: "200"
        content:
          application/json;charset=UTF-8:
            schema:
              $ref: '#/components/schemas/rest-v1-enums-areaEnums-
455423783'
            examples:
              enum/areaEnum:
                value: |-
                  {

```

```

        "AREA_EAST" : "Ida",
        "AREA_WEST" : "Lääne",
        "AREA_SOUTH" : "Lõuna",
        "AREA_NORTH" : "Põhja"
    }
}
/rest/v1/enums/conditionEnums:
get:
  tags:
  - enums
  summary: Kõikide ConditionEnum enumite päring
  description: Kõikide ConditionEnum enumite päring
  operationId: enum/conditionEnum
  parameters:
  - name: Accept
    in: header
    description: Peab olema `application/json`
    required: true
    schema:
      type: string
      default: application/json
      example: application/json
  responses:
    "200":
      description: "200"
      content:
        application/json;charset=UTF-8:
          schema:
            $ref: '#/components/schemas/rest-v1-enums-
conditionEnums1320646597'
          examples:
            enum/conditionEnum:
              value: |-
                {
                  "CONDITION_AND" : "ja",
                  "CONDITION_OR" : "või"
                }
}
/rest/v1/enums/deviceTypeEnum:
get:
  tags:
  - enums
  summary: Kõikide DeviceTypeEnum enumite päring
  description: Kõikide DeviceTypeEnum enumite päring
  operationId: enum/deviceTypeEnum
  parameters:
  - name: Accept
    in: header
    description: Peab olema `application/json`
    required: true
    schema:
      type: string
      default: application/json
      example: application/json
  responses:
    "200":
      description: "200"
      content:
        application/json;charset=UTF-8:
          schema:

```

```

    $ref: '#/components/schemas/rest-v1-enums-
deviceTypeEnums1446031985'
  examples:
    enum/deviceTypeEnum:
      value: |-
        {
          "DEVICE_TYPE_EMAIL" : "E-mail",
          "DEVICE_TYPE_EMAIL_CHIEF" : "E-mail autovanem",
          "DEVICE_TYPE_MGIS" : "MGIS",
          "DEVICE_TYPE_MOBILE" : "MOBIIL",
          "DEVICE_TYPE_MOBILE_CHIEF" : "MOBIIL autovanem",
          "DEVICE_TYPE_NUGIS" : "NUGIS",
          "DEVICE_TYPE_TETRA" : "TETRA",
          "DEVICE_TYPE_TETRA_VEHICLE" : "TETRA autojaam",
          "DEVICE_TYPE_TETRA_CHIEF" : "TETRA autovanem",
          "DEVICE_TYPE_TETRA_PERSON" : "TETRA käsijaam",
          "DEVICE_TYPE_DERIVED" : "Tuletatud"
        }
/rest/v1/enums/eventStatusEnums:
get:
  tags:
  - enums
  summary: Kõikide EventStatusEnum enumite päring
  description: Kõikide EventStatusEnum enumite päring
  operationId: enum/eventStatusEnum
  parameters:
  - name: Accept
    in: header
    description: Peab olema `application/json`
    required: true
    schema:
      type: string
      default: application/json
    example: application/json
  responses:
    "200":
      description: "200"
      content:
        application/json;charset=UTF-8:
          schema:
            $ref: '#/components/schemas/rest-v1-enums-
eventStatusEnums-1396277672'
          examples:
            enum/eventStatusEnum:
              value: |-
                {
                  "STATUS_DELETED" : "Kustutatud",
                  "STATUS_DRAFT" : "Mustand",
                  "STATUS_OPEN" : "Töös",
                  "STATUS_OPEN_WITHOUT_RESOURCES" : "Töös
ressursita",
                  "STATUS_OPEN_PARTLY" : "Töös osaliselt",
                  "STATUS_OPEN_WITH_RESOURCES" : "Töös ressursiga",
                  "STATUS_OPEN_RESOURCES_LEFT" : "Töös ressursid
lahkunud",
                  "STATUS_CLOSED" : "Lõpetatud",
                  "STATUS_ANNULLED" : "Annulleeritud"
                }

```

```

/rest/v1/enums/priorityEnums:
  get:
    tags:
      - enums
    summary: Kõikide PriorityEnums enumite päring
    description: Kõikide PriorityEnums enumite päring
    operationId: enum/priorityEnum
    parameters:
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
          example: application/json
    responses:
      "200":
        description: "200"
        content:
          application/json;charset=UTF-8:
            schema:
              $ref: '#/components/schemas/rest-v1-enums-priorityEnums-
873454254'
        examples:
          enum/priorityEnum:
            value: |-
              {
                "PRIORITY_1" : {
                  "value": "1",
                  "shortValue": "1",
                  "orderNo": 100,
                  "weight": 100,
                  "department": [ "SERVICE_RESCUE" ]
                },
                "PRIORITY_1A" : {
                  "value": "1A",
                  "shortValue": "1A",
                  "orderNo": 90,
                  "weight": 90,
                  "department": [ "SERVICE_RESCUE" ]
                },
                "PRIORITY_2" : {
                  "value": "2",
                  "shortValue": "2",
                  "orderNo": 110,
                  "weight": 110,
                  "department": [ "SERVICE_RESCUE" ]
                },
                "PRIORITY_3" : {
                  "value": "3",
                  "shortValue": "3",
                  "orderNo": 120,
                  "weight": 120,
                  "department": [ "SERVICE_RESCUE" ]
                },
                "PRIORITY_4" : {
                  "value": "4",

```



```

        "shortValue" : "4",
        "orderNo" : 130,
        "weight" : 130,
        "department" : [ "SERVICE_RESCUE" ]
    },
    "PRIORITY_A" : {
        "value" : "A",
        "shortValue" : "A",
        "orderNo" : 140,
        "weight" : 140,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_B" : {
        "value" : "B",
        "shortValue" : "B",
        "orderNo" : 150,
        "weight" : 150,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_C" : {
        "value" : "C",
        "shortValue" : "C",
        "orderNo" : 160,
        "weight" : 160,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE" ]
    },
    "PRIORITY_D" : {
        "value" : "D",
        "shortValue" : "D",
        "orderNo" : 170,
        "weight" : 170,
        "department" : [ "SERVICE_AMBULANCE" ]
    },
    "PRIORITY_INFO" : {
        "value" : "Info",
        "shortValue" : "I",
        "orderNo" : 20,
        "weight" : 20,
        "department" : [ "SERVICE_POLICE" ]
    },
    "PRIORITY_CONSULTATION" : {
        "value" : "Nõuanne",
        "shortValue" : "N",
        "orderNo" : 10,
        "weight" : 10,
        "department" : [ "SERVICE_AMBULANCE",
"SERVICE_POLICE", "SERVICE_RESCUE" ]
    }
}

/rest/v1/enums/resourceStatusEnums:
get:
tags:
- enums
summary: Kõikide ResourceStatusEnum enumite päring
description: Kõikide ResourceStatusEnum enumite päring

```

```

operationId: enum/resourceStatusEnum
parameters:
- name: Accept
  in: header
  description: Peab olema `application/json`
  required: true
  schema:
    type: string
    default: application/json
    example: application/json
responses:
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/rest-v1-enums-
resourceStatusEnums-1756525789'
        examples:
          enum/resourceStatusEnum:
            value: |-
              {
                "RESOURCE_STATUS_ANNULLED" : "Annulleeritud",
                "RESOURCE_STATUS_AWAY_BUSY" : "Eemal hõivatud",
                "RESOURCE_STATUS_AWAY" : "Eemal vaba",
                "RESOURCE_STATUS_AT_HOSPITAL" : "Haiglas",
                "RESOURCE_STATUS_TO_HOSPITAL" : "Haiglasse",
                "RESOURCE_STATUS_AT_HOME" : "Kodus",
                "RESOURCE_STATUS_ON_SITE" : "Kohal",
                "RESOURCE_STATUS_LEAVING_BUSY" : "Lahkun hõivatud",
                "RESOURCE_STATUS_LEAVING" : "Lahkun vaba",
                "RESOURCE_STATUS_ALARMED" : "Teatamine",
                "RESOURCE_STATUS_ALARMED_BY_RADIO" : "Teatamine
raadio teel",
                "RESOURCE_STATUS_ACCEPTED" : "Vastuvõetud",
                "RESOURCE_STATUS_DEPARTURE" : "Väljasõit",
                "RESOURCE_STATUS_DEPARTURE_BY_RADIO" : "Väljasõit
raadio teel"
              }
/rest/v1/resource/alarm:
  post:
    tags:
    - resource
    summary: Ressurside alarmeerimine. Piiratud autenditud kasutajale
(smitTokenId'ga).
    description: Ressurside alarmeerimine. Piiratud autenditud
kasutajale (smitTokenId'ga).
    operationId: alarm/index
    parameters:
    - name: Accept
      in: header
      description: Peab olema `application/json`
      required: true
      schema:
        type: string
        default: application/json
        example: application/json
    - name: smitTokenID

```

```

    in: header
    description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId` 'd
  required: true
  schema:
    type: string
    example: 29c190e8-5113-4cb3-843f-1d4cfa954842
requestBody:
  content:
    text/plain;charset=ISO-8859-1:
      schema:
        $ref: '#/components/schemas/AlarmRequest'
      examples:
        alarm/index-forbidden:
          value: |-
            {
              "alarmedResourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
              "eventNo" : "EVENT_NO",
              "alarmDTime" : "2024-01-30T13:28:08Z",
              "forced" : false
            }
        alarm/index:
          value: |-
            {
              "alarmedResourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
              "eventNo" : "EVENT_NO",
              "alarmDTime" : "2024-01-30T13:28:08Z",
              "forced" : false
            }
        application/json:
          schema:
            $ref: '#/components/schemas/AlarmRequest'
          examples:
            alarm/index-notValid:
              value: |-
                {
                  "alarmedResourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
                  "eventNo" : null,
                  "alarmDTime" : "2024-01-30T13:28:08Z",
                  "forced" : false
                }
      responses:
        "401":
          description: "401"
          content:
            application/json;charset=UTF-8:
              schema:
                $ref: '#/components/schemas/CommonResponse'
              examples:
                alarm/index-forbidden:
                  value: |-
                    {
                      "status" : 401,
                      "message" : "Kasutajal puudub kehtiv sessioon"
                    }

```

```

    }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/AlarmResponse'
        examples:
          alarm/index:
            value: |-
              {
                "alarmedResources" : [ {
                  "name" : "RESOURCE_NAME_1",
                  "status" : "RESOURCE_STATUS_ALARMED",
                  "estimatedOnSiteArrivalTimeInSeconds" : 1282.8,
25402.8
                  "estimatedOnSiteArrivalDistanceInMeters" :
                    } ],
                "notAlarmedResources" : [ {
                  "name" : "RESOURCE_NAME_2",
                  "reason" : "Kutsungile RESOURCE_NAME_2 ei saanud
korraldust anda."
                    } ]
              }
  "412":
    description: "412"
    content:
      application/json;charset=UTF-8:
        schema:
          $ref: '#/components/schemas/CommonResponse'
        examples:
          alarm/index-notValid:
            value: |-
              {
                "status" : 412,
                "message" : "\"eventNo\" peab olema määratud",
                "errorHeader" : "Oleku muutmine ebaõnnestus!"
              }
  /rest/v1/resource/setStatus:
    post:
      tags:
        - resource
      summary: Ressurside staatuse muutmine. Piiratud autenditud
kasutajale (smitTokenId'ga).
      description: Ressurside staatuse muutmine. Piiratud autenditud
kasutajale (smitTokenId'ga).
      operationId: status/index
      parameters:
        - name: Accept
          in: header
          description: Peab olema `application/json`
          required: true
          schema:
            type: string
            default: application/json
          example: application/json
        - name: smitTokenID
          in: header

```

```

description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId` 'd
required: true
schema:
  type: string
example: a0e3037b-a7c4-485c-9c51-625a9428c6cc
requestBody:
content:
  text/plain;charset=ISO-8859-1:
  schema:
    $ref: '#/components/schemas/StatusRequest'
  examples:
    status/index-forbidden:
      value: |-
        {
          "resourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
          "newStatusEnum" : "RESOURCE_STATUS_ANNULLED",
          "statedTime" : "2024-01-30T13:28:08Z",
          "additionalInfo" : "ADDITIONAL_INFO"
        }
    status/index:
      value: |-
        {
          "resourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
          "newStatusEnum" : "RESOURCE_STATUS_ANNULLED",
          "statedTime" : "2024-01-30T13:28:08Z",
          "additionalInfo" : "ADDITIONAL_INFO"
        }
  application/json:
  schema:
    $ref: '#/components/schemas/StatusRequest'
  examples:
    status/index-notValid:
      value: |-
        {
          "resourceNames" : [ "RESOURCE_NAME_1",
"RESOURCE_NAME_2" ],
          "newStatusEnum" : "unknown",
          "statedTime" : "2024-01-30T13:28:08Z",
          "additionalInfo" : "ADDITIONAL_INFO"
        }
responses:
  "401":
description: "401"
content:
  application/json;charset=UTF-8:
  schema:
    $ref: '#/components/schemas/CommonResponse'
  examples:
    status/index-forbidden:
      value: |-
        {
          "status" : 401,
          "message" : "Kasutajal puudub kehtiv sessioon"
        }

```

```

"200":
  description: "200"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/StatusResponse'
      examples:
        status/index:
          value: |-
            {
              "changedResources" : [ {
                "name" : "RESOURCE_NAME_1",
                "status" : "RESOURCE_STATUS_ANNULLED"
              } ],
              "notChangedResources" : [ {
                "name" : "RESOURCE_NAME_2",
                "reason" : "Ressursi RESOURCE_NAME_2 staatust
RESOURCE_STATUS_AT_HOME ei saa muuta annulleerituks!"
              } ]
            }

"412":
  description: "412"
  content:
    application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/CommonResponse'
      examples:
        status/index-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "\"newStatusEnum\" peab olema
määratud",
              "errorHeader" : "Oleku muutmine ebaõnnestus!"
            }

/rest/v1/userSetting/save:
  get:
    tags:
      - userSetting
    summary: "Kasutaja seadete salvestamine. Vajab autenditud kasutajat
(smitTokenId),\
  \ kelle alla seadet salvestada"
    description: "Kasutaja seadete salvestamine. Vajab autenditud
kasutajat (smitTokenId),\
  \ kelle alla seadet salvestada"
    operationId: userSetting/save-forbidden
    parameters:
      - name: typeEnum
        in: query
        description: Seadetele antud tunnus
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:

```

```

    type: string
    default: application/json
    example: application/json
  requestBody:
    content:
      text/plain;charset=ISO-8859-1:
        schema:
          $ref: '#/components/schemas/UserSettingSaveRequest'
        examples:
          userSetting/save-forbidden:
            value: |-
              {
                "filter1" : "45a30158-06fe-47cb-b8bf-5b98991756d6",
                "filter2" : [ "9fd1a010-7a99-4aac-ace2-436e0927b56f",
"6d50472c-b939-453e-88c6-8e5ee475b759", "c331c643-eb8f-4742-a834-
c40d482d12fd" ]
              }
    responses:
      "401":
        description: "401"
        content:
          application/json;charset=UTF-8:
            schema:
              $ref: '#/components/schemas/CommonResponse'
            examples:
              userSetting/save-forbidden:
                value: |-
                  {
                    "status" : 401,
                    "message" : "Kasutajal puudub kehtiv sessioon"
                  }
  post:
    tags:
      - userSetting
    summary: "Kasutaja seadete salvestamine. Vajab autenditud kasutajat (smitTokenId), \
      \ kelle alla seadete salvestada"
    description: "Kasutaja seadete salvestamine. Vajab autenditud kasutajat (smitTokenId), \
      \ kelle alla seadete salvestada"
    operationId: userSetting/save
    parameters:
      - name: typeEnum
        in: query
        description: Seadetele antud tunnus
        required: true
        schema:
          type: string
      - name: Accept
        in: header
        description: Peab olema `application/json`
        required: true
        schema:
          type: string
          default: application/json
          example: application/json
      - name: smitTokenID
        in: header

```

```

description: Antud teenus vajab autentimiseks/autoriseerimiseks
HK portaali
  `smitTokenId` 'd
required: true
schema:
  type: string
example: f37fb5df-4658-48e3-b980-e49a8daf862a
requestBody:
content:
  text/plain;charset=ISO-8859-1:
  schema:
    $ref: '#/components/schemas/rest-v1-userSetting486549215'
  examples:
    userSetting/save-notValid:
      value: |-
        {
          "filter1" : "595dab8f-bdc0-4f37-9fbf-8d1c5da654db",
          "filter2" : [ "f24638cf-eaee-4692-b7e0-7e56ee720ab5",
"8e8348cc-6214-458b-bee8-a3e1c3c18b6a", "efc1c7c2-a660-4625-8979-
3bc1b23eff0d" ]
        }
    userSetting/save:
      value: |-
        {
          "filter1" : "0cc78dc2-58e4-466d-84d6-73d04b9fbaea",
          "filter2" : [ "d7b5dba9-877f-4216-a29e-0b8c4928ac4e",
"f526570f-f5b8-4236-bace-7187b4e2c6bd", "5ab37414-852f-46c1-ac9e-
f6fea639c19e" ]
        }
responses:
  "412":
    description: "412"
    content:
      application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/CommonResponse'
      examples:
        userSetting/save-notValid:
          value: |-
            {
              "status" : 412,
              "message" : "TypeEnum pole lubatud väärtusega,
null"
            }
  "200":
    description: "200"
    content:
      application/json;charset=UTF-8:
      schema:
        $ref: '#/components/schemas/rest-v1-userSetting486549215'
      examples:
        userSetting/save:
          value: |-
            {
              "status" : 202,
              "message" : "OK",
              "userSettings" : {

```



```

        "filter1" : "0cc78dc2-58e4-466d-84d6-
73d04b9fbaea",
        "filter2" : [ "d7b5dba9-877f-4216-a29e-
0b8c4928ac4e", "f526570f-f5b8-4236-bace-7187b4e2c6bd", "5ab37414-852f-
46c1-ac9e-f6fea639c19e" ]
    }
}
components:
  schemas:
    DispatchPlanResponse:
      title: DispatchPlanResponse
      type: object
      properties:
        resourceNeeds:
          type: array
          description: Väljasõiduplaani vajadused
          items:
            type: object
            properties:
              alarmed:
                type: boolean
                description: Vajadusega seotud ressursile on antud
korraldus
              resource:
                type: object
                properties:
                  distanceMeters:
                    type: number
                    description: Eelduslik kohale jõudmise teekonna
pikkus
                  vehicleTypeAmbulanceEnumText:
                    type: string
                    description: Ressursi tehnikaliigi nimetus
                  lastTimeOnEvent:
                    type: string
                    description: Viimati sündmusel oldud aeg ISO
formaadis - `yyyy-MM-dd'T'HH:mm:ssXX`
                  vehicleTypeAmbulanceEnum:
                    type: string
                    description: Ressursi tehnikaliigi kood.
VehicleTypeAmbulanceEnum
                    väärtused on klassikaator tüüpi koodid Reegel
rakendusest
                    name:
                      type: string
                      description: Ressursi kutsung
                    externalId:
                      type: string
                      description: "Ressursi väline unikaalne tunnus, ID"
                    distanceSeconds:
                      type: number
                      description: Eelduslik kohale jõudmise aeg sekundites
                    id:
                      type: number
                      description: "Ressursi identifikaator rakenduses,
ID"
                    stationId:
                      type: string

```

```

        description: "Ressursi komando/üksuse tunnus, ID"
        description: Vajaduse katmiseks leitud ressurss
cancelled:
    type: boolean
    description: Vajadus on märgitud mitte vajalikuks
active:
    type: boolean
    description: Vajadus on aktiiven
extraResourceVehicleTypeEnum:
    type: string
    description: Lisa vajaduse tehnikaliigi kood. Väärtused
on klassikaator
    tüüpi koodid Reegel rakendusest
done:
    type: boolean
    description: Vajadus on täidetud (seotud ressurss kohal
käinud)
ruleResourceNeed:
    type: object
    properties:
        ruleResourceNeedTerms:
            type: array
            description: Reegli tingimused
            items:
                type: object
                properties:
                    orderNo:
                        type: number
                        description: Reegli tingimuse jrk. nr
                    terms:
                        type: object
                        properties:
                            vehicleTypeCode:
                                type: string
                                description: ""
                            description: Reegli tingimuse kirjeldus Reegel
rakendusest
amount:
    type: number
    description: Reeglist tulenev vajaduse kogus
orderNo:
    type: number
    description: Reegli jrk. nr
conditionEnum:
    type: string
    description: "Reegli tingimuste tüübi kood.
ConditionEnum vää\
        rtused - CONDITION_AND, CONDITION_OR"
        description: Vajadus tuleneb reeglist
extraResourceVehicleTypeEnumText:
    type: string
    description: Lisa vajaduse tehnikaliigi nimetus
eventVersion:
    type: number
    description: "Väljasõiduplaaniga seotud sündmuse versioon,
millelt plaan\
        \ on arvutatud"
dispatchRuleName:

```

```

    type: string
    description: Väljasõiduplaani koostamise aluseks oleva reegli
nimetus
    preferredResourceStationIds:
      type: string
      description: Väljasõiduplaani reeglist tulenevad eelistatud
üksuste/komando
        tunnuste komadega eraldatud tunnuste nimekiri
    dispatchRuleId:
      type: string
      description: "Väljasõiduplaani koostamise aluseks oleva reegli
tunnus, ID"
    alternativeResources:
      type: array
      description: Väljasõiduplaani alternatiivsed ressursid
      items:
        type: object
        properties:
          distanceMeters:
            type: number
            description: Eelduslik kohale jõudmise teekonna pikkus
          vehicleTypeAmbulanceEnumText:
            type: string
            description: Ressursi tehnikaliigi nimetus
          lastTimeOnEvent:
            type: string
            description: Viimati sündmusel oldud aeg ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
          vehicleTypeAmbulanceEnum:
            type: string
            description: Ressursi tehnikaliigi kood.
VehicleTypeAmbulanceEnum
          väärtused on klassikaator tüüpi koodid Reegel
rakendusest
            name:
              type: string
              description: Ressursi kutsung
            externalId:
              type: string
              description: "Ressursi väline unikaalne tunnus, ID"
            distanceSeconds:
              type: number
              description: Eelduslik kohale jõudmise aeg sekundites
            id:
              type: number
              description: "Ressursi identifikaator rakenduses, ID"
            stationId:
              type: string
              description: "Ressursi komando/üksuse tunnus, ID"
          eventNo:
            type: string
            description: Väljasõiduplaaniga seotud sündmuse number
          dispatchPlanId:
            type: number
            description: "Väljasõiduplaani identifikaator, ID"
rest-v1-enums-deviceTypeEnums1446031985:
  type: object
  properties:

```

```
DEVICE_TYPE_TETRA_PERSON:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_MGIS:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_NUGIS:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_MOBILE_CHIEF:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_MOBILE:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_EMAIL:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_EMAIL_CHIEF:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_DERIVED:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_TETRA_VEHICLE:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_TETRA:
  type: string
  description: Enum'i nimi ja väärtus
DEVICE_TYPE_TETRA_CHIEF:
  type: string
  description: Enum'i nimi ja väärtus
rest-v1-enums-priorityEnums-873454254:
  type: object
  properties:
    PRIORITY_C:
      type: object
      properties:
        orderNo:
          type: number
          description: Järjekorra number nimekirjades
        weight:
          type: number
          description: Prioriteedi kaal suhtes teiste prioriteetidega
        shortValue:
          type: string
          description: Prioriteedi lühinimi
        department:
          type: array
          description: "Teenus(ed), kus kasutuses"
          items:
            oneOf:
              - type: object
              - type: boolean
              - type: string
              - type: number
        value:
```

```
    type: string
    description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_3:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste prioriteetidega
    shortValue:
      type: string
      description: Prioriteedi lühinimi
    department:
      type: array
      description: "Teenus(ed), kus kasutuses"
      items:
        oneOf:
          - type: object
          - type: boolean
          - type: string
          - type: number
    value:
      type: string
      description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_D:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste prioriteetidega
    shortValue:
      type: string
      description: Prioriteedi lühinimi
    department:
      type: array
      description: "Teenus(ed), kus kasutuses"
      items:
        oneOf:
          - type: object
          - type: boolean
          - type: string
          - type: number
    value:
      type: string
      description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_4:
  type: object
  properties:
    orderNo:
      type: number
```

```
    description: Järjekorra number nimekirjades
weight:
  type: number
  description: Prioriteedi kaal suhtes teiste prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
  type: array
  description: "Teenus(ed), kus kasutuses"
  items:
    oneOf:
      - type: object
      - type: boolean
      - type: string
      - type: number
value:
  type: string
  description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_1A:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
weight:
  type: number
  description: Prioriteedi kaal suhtes teiste prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
  type: array
  description: "Teenus(ed), kus kasutuses"
  items:
    oneOf:
      - type: object
      - type: boolean
      - type: string
      - type: number
value:
  type: string
  description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_INFO:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
weight:
  type: number
  description: Prioriteedi kaal suhtes teiste prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
```

```
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
  value:
    type: string
    description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_A:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste prioriteetidega
    shortValue:
      type: string
      description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
  value:
    type: string
    description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_1:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste prioriteetidega
    shortValue:
      type: string
      description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
```

```
value:
  type: string
  description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_B:
  type: object
properties:
  orderNo:
    type: number
    description: Järjekorra number nimekirjades
  weight:
    type: number
    description: Prioriteedi kaal suhtes teiste prioriteetidega
  shortValue:
    type: string
    description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
value:
  type: string
  description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_2:
  type: object
properties:
  orderNo:
    type: number
    description: Järjekorra number nimekirjades
  weight:
    type: number
    description: Prioriteedi kaal suhtes teiste prioriteetidega
  shortValue:
    type: string
    description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
value:
  type: string
  description: Prioriteedi nimetus
description: Enum'i nimi ja väärtus
PRIORITY_CONSULTATION:
  type: object
properties:
  orderNo:
```



```
    type: number
    description: Järjekorra number nimekirjades
weight:
  type: number
  description: Prioriteedi kaal suhtes teiste prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
  type: array
  description: "Teenus(ed), kus kasutuses"
  items:
    oneOf:
      - type: object
      - type: boolean
      - type: string
      - type: number
value:
  type: string
  description: Prioriteedi nimetus
  description: Enum'i nimi ja väärtus
rest-v1-enums-areaEnums-455423783:
  type: object
  properties:
    AREA_SOUTH:
      type: string
      description: Enum'i nimi ja väärtus
    AREA_NORTH:
      type: string
      description: Enum'i nimi ja väärtus
    AREA_EAST:
      type: string
      description: Enum'i nimi ja väärtus
    AREA_WEST:
      type: string
      description: Enum'i nimi ja väärtus
AlarmRequest:
  title: AlarmRequest
  type: object
  properties:
    alarmDTime:
      type: string
      description: Alameerimise aeg ISO formaadis - `yyyy-MM-dd'T'HH:mm:ssXX`
    alarmedResourceNames:
      type: array
      description: Ressurside kutsungid
      items:
        oneOf:
          - type: object
          - type: boolean
          - type: string
          - type: number
forced:
  type: boolean
  description: Kas sundida alameerimine
eventNo:
  type: string
```

```

        description: Sündmuse number
HospitalResponse:
  title: HospitalResponse
  type: object
  properties:
    total:
      type: number
      description: Kogu kirjete arv. Tagastatud kirjete arv võib olla
sellest
      väiksem - tuleb pärida lehekülgede kaupa
  list:
    type: array
    description: "Päringu tulemus, nimekiri päritud objektidest.
Haigla mudel"
  items:
    type: object
    properties:
      address:
        type: string
        description: Haigla aadress
      code:
        type: string
        description: Haigla kood
      areaEnumText:
        type: string
        description: Piirkonna nimetus
      areaEnum:
        type: string
        description: "Haigla piirkonna kood (HK vaates). AreaEnum
väärtused\
          \ - AREA_NORTH, AREA_SOUTH, AREA_EAST, AREA_WEST"
      description:
        type: string
        description: Lisa kirjeldus
      editorSessionName:
        type: string
        description: Viimase muutja täisnimi
      creatorSessionName:
        type: string
        description: Looja täisnimi
      addressWgsLongitude:
        type: number
        description: Asukoha pikkuskraad WGS84 projektsioonis
      version:
        type: number
        description: Haigla kirje versioon
      lastUpdated:
        type: string
        description: Muutmise kuupäeva ISO formaadis - `yyyy-MM-
dd'T'HH:mm:ssXX`
      editorSessionIdCode:
        type: string
        description: Viimase muutja isikukood
      creatorSessionIdCode:
        type: string
        description: Looja isikukood
      dateCreated:
        type: string

```

```
description: Loomise kuupäev ISO formaadis - `yyyy-MM-  
dd'T'HH:mm:ssXX`  
deleted:  
  type: boolean  
  description: Kas on kustutatud  
name:  
  type: string  
  description: Haigla nimetus  
id:  
  type: number  
  description: "Haigla identifikaator rakenduses, ID"  
addressWgsLatitude:  
  type: number  
  description: Asukoha laiuskraad WGS84 projektsioonis  
CommonResponse:  
  title: CommonResponse  
  type: object  
  properties:  
    instruction:  
      type: string  
      description: "Kui olemas, siis juhised lõppkasutajale, kuidas  
jätkata"  
    errorCode:  
      type: string  
      description: "Kui olemas, siis veakood"  
    cause:  
      type: string  
      description: "Kui olemas, siis põhjus originaal veast,  
tehniline teade"  
    message:  
      type: string  
      description: "OK, kui päring õnnestus või ebaõnnestumise vea  
teade"  
    errorHeader:  
      type: string  
      description: "Kui olemas, siis veateate pealdis"  
    status:  
      type: number  
      description: "HTTP vastuse oleku kood (status code), 202 -  
succeeded"  
rest-v1-enums-resourceStatusEnums-1756525789:  
  type: object  
  properties:  
    RESOURCE_STATUS_ANNULLED:  
      type: string  
      description: Enum'i nimi ja väärtus  
    RESOURCE_STATUS_LEAVING:  
      type: string  
      description: Enum'i nimi ja väärtus  
    RESOURCE_STATUS_AT_HOME:  
      type: string  
      description: Enum'i nimi ja väärtus  
    RESOURCE_STATUS_AWAY:  
      type: string  
      description: Enum'i nimi ja väärtus  
    RESOURCE_STATUS_TO_HOSPITAL:  
      type: string  
      description: Enum'i nimi ja väärtus
```

```

RESOURCE_STATUS_ACCEPTED:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_ON_SITE:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_DEPARTURE_BY_RADIO:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_DEPARTURE:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_AWAY_BUSY:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_ALARMED:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_LEAVING_BUSY:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_ALARMED_BY_RADIO:
  type: string
  description: Enum'i nimi ja väärtus
RESOURCE_STATUS_AT_HOSPITAL:
  type: string
  description: Enum'i nimi ja väärtus
StatusResponse:
  title: StatusResponse
  type: object
  properties:
    changedResources:
      type: array
      description: Muudetud ressursid
      items:
        type: object
        properties:
          name:
            type: string
            description: Muudetud ressursi nimi
          status:
            type: string
            description: "Muudetud ressursi staatus, <<enums-
resourceStatusEnum>>\
          \ väärtus - RESOURCE_STATUS_ALARMED,
RESOURCE_STATUS_ALARMED_BY_RADIO,\
          \ RESOURCE_STATUS_ON_SITE, RESOURCE_STATUS_AT_HOME,
RESOURCE_STATUS_DEPARTURE,\
          \ RESOURCE_STATUS_DEPARTURE_BY_RADIO,
RESOURCE_STATUS_TO_HOSPITAL,\
          \ RESOURCE_STATUS_AT_HOSPITAL, RESOURCE_STATUS_AWAY,
RESOURCE_STATUS_AWAY_BUSY,\
          \ RESOURCE_STATUS_LEAVING,
RESOURCE_STATUS_LEAVING_BUSY, RESOURCE_STATUS_ANNULLED,\
          \ RESOURCE_STATUS_ACCEPTED"
    notChangedResources:
      type: array
      description: Mittemuudetud ressursid

```

```

    items:
      type: object
      properties:
        reason:
          type: string
          description: Mittemuudetud ressursi põhjus
        name:
          type: string
          description: Mittemuudetud ressursi nimi
StatusRequest:
  title: StatusRequest
  type: object
  properties:
    resourceNames:
      type: array
      description: Ressurside kutsungid
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
    stateDTime:
      type: string
      description: Ressursi oleku aeg ISO formaadis - `yyyy-MM-dd'T'HH:mm:ssXX`
    newStatusEnum:
      type: string
      description: "Ressursi uus olek, <<enums-resourceStatusEnum>>
väärtus -\
  \ RESOURCE_STATUS_ALARMED, RESOURCE_STATUS_ALARMED_BY_RADIO,
RESOURCE_STATUS_ON_SITE,\
  \ RESOURCE_STATUS_AT_HOME, RESOURCE_STATUS_DEPARTURE,
RESOURCE_STATUS_DEPARTURE_BY_RADIO,\
  \ RESOURCE_STATUS_TO_HOSPITAL, RESOURCE_STATUS_AT_HOSPITAL,
RESOURCE_STATUS_AWAY,\
  \ RESOURCE_STATUS_AWAY_BUSY, RESOURCE_STATUS_LEAVING,
RESOURCE_STATUS_LEAVING_BUSY,\
  \ RESOURCE_STATUS_ANNULLED, RESOURCE_STATUS_ACCEPTED"
    additionalInfo:
      type: string
      description: Lisainfo
    hospitalCode:
      type: string
      description: Ressursi uue haigla kood
    eventNo:
      type: string
      description: Oleku muudatusega seotud sündmuse number
rest-v1-enums-conditionEnums1320646597:
  type: object
  properties:
    CONDITION_AND:
      type: string
      description: Enum'i nimi ja väärtus
    CONDITION_OR:
      type: string
      description: Enum'i nimi ja väärtus
ResourceResponse:

```

```

title: ResourceResponse
type: object
properties:
  total:
    type: number
    description: Kogu kirjete arv. Tagastatud kirjete arv võib olla
sellest
    väiksem - tuleb pärida lehekülgede kaupa
  list:
    type: array
    description: "Päringu tulemus, nimekiri päritud objektidest.
Ressursi mudel"
    items:
      type: object
      properties:
        expiring:
          type: boolean
          description: Kas ressurss on aegumas
        areaEnumText:
          type: string
          description: Ressursi piirkonna nimetus
        areaEnum:
          type: string
          description: "Ressursi piirkonna kood. AreaEnum väärtused
- AREA_NORTH, \
    \ AREA_SOUTH, AREA_EAST, AREA_WEST"
        vehicleTypeAmbulanceEnum:
          type: string
          description: Ressursi tehnikaliigi kood.
VehicleTypeAmbulanceEnum
          väärtused on klassikaator tüüpi koodid Reegel
rakendusest
        stationLatitude:
          type: number
          description: Ressursi komando/üksuse laiuskraad WGS84
projektsioonis
        stationLongitude:
          type: number
          description: Ressursi komando/üksuse pikkuskraad WGS84
projektsioonis
        creatorSessionName:
          type: string
          description: Looja täisnimi
        responseRange:
          type: number
          description: Reageerimise kaugus kilomeetrites
        lastUpdated:
          type: string
          description: Muutmise kuupäeva ISO formaadis - `yyyy-MM-
dd'T'HH:mm:ssXX`
        creatorSessionIdCode:
          type: string
          description: Looja isikukood
        dateCreated:
          type: string
          description: Loomise kuupäev ISO formaadis - `yyyy-MM-
dd'T'HH:mm:ssXX`
        members:

```

```

type: array
items:
  type: object
  properties:
    lastUpdated:
      type: string
      description: Muutmise kuupäeva ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
    editorSessionIdCode:
      type: string
      description: Viimase muutja isikukood
    creatorSessionIdCode:
      type: string
      description: Looja isikukood
    resourceId:
      type: number
      description: "Viide ressursile, ID"
    dateCreated:
      type: string
      description: Loomise kuupäev ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
    isChief:
      type: boolean
      description: Kas on meeskonnavanem
    active:
      type: boolean
      description: Kas liige on aktiivne
    externalId:
      type: string
      description: Ressurss liikme väline tunnus
    id:
      type: number
      description: "Ressursi liikme identifikaator
rakenduses, ID"
    editorSessionName:
      type: string
      description: Viimase muutja täisnimi
    creatorSessionName:
      type: string
      description: Looja täisnimi
    version:
      type: number
      description: Ressursi liikme kirje versioon
    station:
      type: string
      description: Ressursi komando/üksuse nimetus
    id:
      type: number
      description: "Resssursi identifikaator rakenduses, ID"
    events:
      type: array
      items:
        type: object
        properties:
          hospitalDate:
            type: string
            description: Haiglasse aeg ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`

```

```

resourceId:
  type: number
  description: "Viide ressursile, ID"
onSiteDate:
  type: string
  description: Kohal aeg ISO formaadis - `yyyy-MM-
dd'T'HH:mm:ssXX`
creatorSessionName:
  type: string
  description: Looja täisnimi
statusEnumText:
  type: string
  description: Ressursi oleku nimetus sündmuse suhtes
acceptedDate:
  type: string
  description: Vastuvõtmise aeg ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
radio:
  type: boolean
  description: Kas alarmeerimine toimus raadio teel
või komandost/üksusest
homeDate:
  type: string
  description: Kodus (komandos/üksuses) aeg ISO
formaadis - `yyyy-MM-dd'T'HH:mm:ssXX`
lastUpdated:
  type: string
  description: Muutmise kuupäeva ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
creatorSessionIdCode:
  type: string
  description: Looja isikukood
statusEnum:
  type: string
  description: "Ressursi oleku kood sündmuse suhtes.
ResourceStatusEnum\
    \ väärtused - RESOURCE_STATUS_ALARMED,
RESOURCE_STATUS_ALARMED_BY_RADIO,\
    \ RESOURCE_STATUS_ON_SITE,
RESOURCE_STATUS_AT_HOME, RESOURCE_STATUS_DEPARTURE,\
    \ RESOURCE_STATUS_DEPARTURE_BY_RADIO,
RESOURCE_STATUS_TO_HOSPITAL,\
    \ RESOURCE_STATUS_AT_HOSPITAL,
RESOURCE_STATUS_AWAY, RESOURCE_STATUS_AWAY_BUSY,\
    \ RESOURCE_STATUS_LEAVING,
RESOURCE_STATUS_LEAVING_BUSY, RESOURCE_STATUS_ANNULLED,\
    \ RESOURCE_STATUS_ACCEPTED"
dateCreated:
  type: string
  description: Loomise kuupäev ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
hospitalId:
  type: number
  description: "Viide haiglale, ID"
additionalData:
  type: string
  description: Lisa andmed/märkused
departureDate:

```



```

        type: string
        description: Väljasõit aeg ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
        id:
            type: number
            description: "Ressursi-sündmuse seose
identifikaator rakenduses, \
        \ ID"
        leavingDate:
            type: string
            description: Lahkumise aeg ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
        stationAlarmedStartDate:
            type: string
            description: Komandos/üksuses alarmeerimise alguse
aeg ISO formaadis
                - `yyyy-MM-dd'T'HH:mm:ssXX`
        eventId:
            type: number
            description: "Viide sündmusele, ID"
        annulledDate:
            type: string
            description: Annulleerimise aeg ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
        alarmedDate:
            type: string
            description: Alarmeerimise aeg ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
        stationAlarmedEndDate:
            type: string
            description: Komandos/üksuses alarmeerimise lõpu
aeg ISO formaadis
                - `yyyy-MM-dd'T'HH:mm:ssXX`
        externalId:
            type: string
            description: "Ressurss-sündmuse seose väline
tunnus, ID"
        active:
            type: boolean
            description: Ressurss-sündmuse seose
aktiivsus/kehtivus
        removedFromSync:
            type: boolean
            description: Seotud andmete alglaadimisega SOSst
        hospitalName:
            type: string
            description: Haigla nimetus
        editorSessionName:
            type: string
            description: Viimase muutja täisnimi
        estimatedOnSiteArrivalTimeInSeconds:
            type: number
            description: Alarmeerimisest eelduslik kohale
jõudmise aeg sekundites
        version:
            type: number
            description: Ressursi-sündmuse seose kirje versioon
        editorSessionIdCode:

```

```

        type: string
        description: Viimase muutja isikukood
stationId:
    type: string
    description: "Ressursi komando/üksuse tunnus, ID"
vehicleTypeAmbulanceEnumText:
    type: string
    description: Ressursi tehnikaliigi nimetus
devices:
    type: array
    items:
        type: object
        properties:
            lastUpdated:
                type: string
                description: Muutmise kuupäeva ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
            resourceId:
                type: number
                description: "Viide seadmega seotud ressursile, ID"
            dateCreated:
                type: string
                description: Loomise kuupäev ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
            externalId:
                type: string
                description: "Seadme väline tunnus - email-i
address, TETRA\
                \ number jne"
            resourceMemberId:
                type: number
                description: "Viide seadmega seotud meeskonna
liikmele, ID"
            id:
                type: number
                description: "Seadme identifikaator rakenduses, ID"
            deviceTypeEnumText:
                type: string
                description: Seadme tüübi nimetus
            version:
                type: number
                description: Seadme kirje versioon
            deviceLocation:
                type: object
                properties:
                    lastUpdated:
                        type: string
                        description: Muutmise kuupäeva ISO formaadis -
`yyyy-MM-dd'T'HH:mm:ssXX`
                    locationTime:
                        type: string
                        description: Asukoha aeg ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
                    latitude:
                        type: number
                        description: Asukoha laiuskraad WGS84
projektsioonis
            id:

```

```

        type: number
        description: "Seadme asukoha identifikaator
rakenduses, \
        \ ID"
    version:
        type: number
        description: Seadme asukoha kirje versioon
    speed:
        type: number
        description: Liikumise kiirus
    direction:
        type: number
        description: Liikumise suund
    longitude:
        type: number
        description: Asukoha pikkuskraad WGS84
projektsioonis
    deviceTypeEnum:
        type: string
        description: "Seadme tüüp. DeviceTypeEnum väärtused
- DEVICE_TYPE_NUGIS, \
        \ DEVICE_TYPE_MGIS, DEVICE_TYPE_TETRA,
DEVICE_TYPE_TETRA_VEHICLE, \
        \ DEVICE_TYPE_TETRA_CHIEF,
DEVICE_TYPE_TETRA_PERSON, DEVICE_TYPE_MOBILE, \
        \ DEVICE_TYPE_MOBILE_CHIEF, DEVICE_TYPE_EMAIL,
DEVICE_TYPE_EMAIL_CHIEF, \
        \ DEVICE_TYPE_DERIVED"
    responseTime:
        type: number
        description: Reageerimise kiirus komandost/üksusest
sekundites
    externalId:
        type: string
        description: "Ressursi väline unikaalne tunnus, ID"
    active:
        type: boolean
        description: Kas ressurss on aktiivne
    editorSessionName:
        type: string
        description: Viimase muutja täisnimi
    resourceLocation:
        type: object
        properties:
            lastUpdated:
                type: string
                description: Muutmise kuupäeva ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
            locationTime:
                type: string
                description: Asukoha aeg ISO formaadis - `yyyy-MM-
dd'T'HH:mm:ssXX`
            latitude:
                type: number
                description: Asukoha laiuskraad WGS84 projektsioonis
    id:
        type: number

```

```

        description: "Ressursi asukoha identifikaator
rakenduses, ID"
        deviceTypeEnumText:
            type: string
            description: Seadme tüüp/asukoha allika nimetus
        version:
            type: number
            description: Ressursi asukoha kirje versioon
        speed:
            type: number
            description: Liikumise kiirus
        direction:
            type: number
            description: Liikumise suund
        longitude:
            type: number
            description: Asukoha pikkuskraad WGS84 projektsioonis
        deviceTypeEnum:
            type: string
            description: "Seadme tüüp/asukoha allikas.
DeviceTypeEnum vää\
            rtused - DEVICE_TYPE_NUGIS, DEVICE_TYPE_MGIS,
DEVICE_TYPE_TETRA,\
            \ DEVICE_TYPE_TETRA_VEHICLE,
DEVICE_TYPE_TETRA_CHIEF, DEVICE_TYPE_TETRA_PERSON,\
            \ DEVICE_TYPE_DERIVED"
        version:
            type: number
            description: Ressursi kirje versioon
        editorSessionIdCode:
            type: string
            description: Viimase muutja isikukood
        resourceStatus:
            type: object
            properties:
                lastUpdated:
                    type: string
                    description: Muutmise kuupäeva ISO formaadis - `yyyy-
MM-dd'T'HH:mm:ssXX`
                editorSessionIdCode:
                    type: string
                    description: Viimase muutja isikukood
                statusNote:
                    type: string
                    description: Oleku lisa märkus
                lastEventDate:
                    type: string
                    description: Viimati sündmusel oldud aeg ISO
formaadis - `yyyy-MM-dd'T'HH:mm:ssXX`
                resourceStatusEnumText:
                    type: string
                    description: Ressursi oleku nimetus
                resourceStatusEnum:
                    type: string
                    description: "Ressursi oleku kood. ResourceStatusEnum
väärtused\
                    \ - RESOURCE_STATUS_ALARMED,
RESOURCE_STATUS_ALARMED_BY_RADIO,\

```

```

        \ RESOURCE_STATUS_ON_SITE, RESOURCE_STATUS_AT_HOME,
RESOURCE_STATUS_DEPARTURE, \
        \ RESOURCE_STATUS_DEPARTURE_BY_RADIO,
RESOURCE_STATUS_TO_HOSPITAL, \
        \ RESOURCE_STATUS_AT_HOSPITAL,
RESOURCE_STATUS_AWAY, RESOURCE_STATUS_AWAY_BUSY, \
        \ RESOURCE_STATUS_LEAVING,
RESOURCE_STATUS_LEAVING_BUSY, RESOURCE_STATUS_ANNULLED, \
        \ RESOURCE_STATUS_ACCEPTED"
    id:
      type: number
      description: "Ressursi oleku identifikaator
rakenduses, ID"
    editorSessionName:
      type: string
      description: Viimase muutja täisnimi
    version:
      type: number
      description: Ressursi oleku kirje versioon
    name:
      type: string
      description: Ressursi kutsung
    validUntil:
      type: string
      description: "Ressurss on kehtiv kuni, aeg ISO formaadis
- `yyyy-MM-dd'T'HH:mm:ssXX`"
    comment:
      type: string
      description: Kommentaar/märkus
    sosVersion:
      type: number
      description: Ressursi versioon SOS'is
AlarmResponse:
  title: AlarmResponse
  type: object
  properties:
    alarmedResources:
      type: array
      description: Alameeritud ressursid
      items:
        type: object
        properties:
          estimatedOnSiteArrivalDistanceInMeters:
            type: number
            description: Alameeritud ressursi läbitav distant
          name:
            type: string
            description: Alameeritud ressursi nimi
          estimatedOnSiteArrivalTimeInSeconds:
            type: number
            description: Alameeritud ressursi saabumise aeg
sekundites
          status:
            type: string
            description: "Alameeritud ressursi staatus, <<enums-
resourceStatusEnum>>\
        \ väärtus - RESOURCE_STATUS_ALARMED,
RESOURCE_STATUS_ALARMED_BY_RADIO, \

```

```
        \ RESOURCE_STATUS_ON_SITE, RESOURCE_STATUS_AT_HOME,
RESOURCE_STATUS_DEPARTURE, \
        \ RESOURCE_STATUS_DEPARTURE_BY_RADIO,
RESOURCE_STATUS_TO_HOSPITAL, \
        \ RESOURCE_STATUS_AT_HOSPITAL, RESOURCE_STATUS_AWAY,
RESOURCE_STATUS_AWAY_BUSY, \
        \ RESOURCE_STATUS_LEAVING,
RESOURCE_STATUS_LEAVING_BUSY, RESOURCE_STATUS_ANNULLED, \
        \ RESOURCE_STATUS_ACCEPTED"
```

notAlarmedResources:

type: array

description: Mitteamarmeeritud ressursid

items:

type: object

properties:

reason:

type: string

description: Mitteamarmeeritud ressursi põhjus

name:

type: string

description: Mitteamarmeeritud ressursi nimi

UserSettingSaveRequest:

title: UserSettingSaveRequest

type: object

rest-v1-enums-eventStatusEnums-1396277672:

type: object

properties:

STATUS_OPEN:

type: string

description: Enum'i nimi ja väärtus

STATUS_DRAFT:

type: string

description: Enum'i nimi ja väärtus

STATUS_OPEN_RESOURCES_LEFT:

type: string

description: Enum'i nimi ja väärtus

STATUS_OPEN_PARTLY:

type: string

description: Enum'i nimi ja väärtus

STATUS_DELETED:

type: string

description: Enum'i nimi ja väärtus

STATUS_OPEN_WITHOUT_RESOURCES:

type: string

description: Enum'i nimi ja väärtus

STATUS_ANNULLED:

type: string

description: Enum'i nimi ja väärtus

STATUS_CLOSED:

type: string

description: Enum'i nimi ja väärtus

STATUS_OPEN_WITH_RESOURCES:

type: string

description: Enum'i nimi ja väärtus

rest-v1-userSetting486549215:

type: object

rest-v1-enums1492299693:

type: object

```

properties:
  resourceStatusEnums:
    type: object
    properties:
      RESOURCE_STATUS_ANNULLED:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_ANNULLED'
      RESOURCE_STATUS_LEAVING:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_LEAVING'
      RESOURCE_STATUS_AT_HOME:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_AT_HOME'
      RESOURCE_STATUS_AWAY:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_AWAY'
      RESOURCE_STATUS_TO_HOSPITAL:
        type: string
        description: ResourceStatusEnum
'RESOURCE_STATUS_TO_HOSPITAL'
      RESOURCE_STATUS_ACCEPTED:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_ACCEPTED'
      RESOURCE_STATUS_ON_SITE:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_ON_SITE'
      RESOURCE_STATUS_DEPARTURE_BY_RADIO:
        type: string
        description: ResourceStatusEnum
'RESOURCE_STATUS_DEPARTURE_BY_RADIO'
      RESOURCE_STATUS_DEPARTURE:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_DEPARTURE'
      RESOURCE_STATUS_AWAY_BUSY:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_AWAY_BUSY'
      RESOURCE_STATUS_ALARMED:
        type: string
        description: ResourceStatusEnum 'RESOURCE_STATUS_ALARMED'
      RESOURCE_STATUS_LEAVING_BUSY:
        type: string
        description: ResourceStatusEnum
'RESOURCE_STATUS_LEAVING_BUSY'
      RESOURCE_STATUS_ALARMED_BY_RADIO:
        type: string
        description: ResourceStatusEnum
'RESOURCE_STATUS_ALARMED_BY_RADIO'
      RESOURCE_STATUS_AT_HOSPITAL:
        type: string
        description: ResourceStatusEnum
'RESOURCE_STATUS_AT_HOSPITAL'
        description: "ResourceStatusEnums - RESOURCE_STATUS_ALARMED,
RESOURCE_STATUS_ALARMED_BY_RADIO, \
  \ RESOURCE_STATUS_ON_SITE, RESOURCE_STATUS_AT_HOME,
RESOURCE_STATUS_DEPARTURE, \
  \ RESOURCE_STATUS_DEPARTURE_BY_RADIO,
RESOURCE_STATUS_TO_HOSPITAL, RESOURCE_STATUS_AT_HOSPITAL, \

```

```

        \ RESOURCE_STATUS_AWAY, RESOURCE_STATUS_AWAY_BUSY,
RESOURCE_STATUS_LEAVING,\
        \ RESOURCE_STATUS_LEAVING_BUSY, RESOURCE_STATUS_ANNULLED,
RESOURCE_STATUS_ACCEPTED"
conditionEnums:
  type: object
  properties:
    CONDITION_AND:
      type: string
      description: ConditionEnum 'CONDITION_AND'
    CONDITION_OR:
      type: string
      description: ConditionEnum 'CONDITION_OR'
  description: "ConditionEnums - CONDITION_AND, CONDITION_OR"
eventStatusEnums:
  type: object
  properties:
    STATUS_OPEN:
      type: string
      description: EventStatusEnum 'STATUS_OPEN'
    STATUS_DRAFT:
      type: string
      description: EventStatusEnum 'STATUS_DRAFT'
    STATUS_OPEN_RESOURCES_LEFT:
      type: string
      description: EventStatusEnum 'STATUS_OPEN_RESOURCES_LEFT'
    STATUS_OPEN_PARTLY:
      type: string
      description: EventStatusEnum 'STATUS_OPEN_PARTLY'
    STATUS_DELETED:
      type: string
      description: EventStatusEnum 'STATUS_DELETED'
    STATUS_OPEN_WITHOUT_RESOURCES:
      type: string
      description: EventStatusEnum
'STATUS_OPEN_WITHOUT_RESOURCES'
    STATUS_ANNULLED:
      type: string
      description: EventStatusEnum 'STATUS_ANNULLED'
    STATUS_CLOSED:
      type: string
      description: EventStatusEnum 'STATUS_CLOSED'
    STATUS_OPEN_WITH_RESOURCES:
      type: string
      description: EventStatusEnum 'STATUS_OPEN_WITH_RESOURCES'
  description: "EventStatusEnums - STATUS_DELETED, STATUS_DRAFT,
STATUS_OPEN,\
        \ STATUS_OPEN_WITHOUT_RESOURCES, STATUS_OPEN_PARTLY,
STATUS_OPEN_WITH_RESOURCES,\
        \ STATUS_OPEN_RESOURCES_LEFT, STATUS_CLOSED, STATUS_ANNULLED"
priorityEnums:
  type: object
  properties:
    PRIORITY_C:
      type: object
      properties:
        orderNo:
          type: number

```


description: Järjekorra number nimekirjades
weight:
type: number
description: Prioriteedi kaal suhtes teiste
prioriteetidega

shortValue:
type: string
description: Prioriteedi lühinimi
department:
type: array
description: "Teenus(ed), kus kasutuses"
items:
oneOf:
- type: object
- type: boolean
- type: string
- type: number

value:
type: string
description: Prioriteedi nimetus
description: PriorityEnum 'PRIORITY_C'

PRIORITY_3:
type: object

properties:
orderNo:
type: number
description: Järjekorra number nimekirjades
weight:
type: number
description: Prioriteedi kaal suhtes teiste

prioriteetidega

shortValue:
type: string
description: Prioriteedi lühinimi
department:
type: array
description: "Teenus(ed), kus kasutuses"
items:
oneOf:
- type: object
- type: boolean
- type: string
- type: number

value:
type: string
description: Prioriteedi nimetus
description: PriorityEnum 'PRIORITY_3'

PRIORITY_D:
type: object

properties:
orderNo:
type: number
description: Järjekorra number nimekirjades
weight:
type: number
description: Prioriteedi kaal suhtes teiste

prioriteetidega

shortValue:

```

    type: string
    description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
  value:
    type: string
    description: Prioriteedi nimetus
  description: PriorityEnum 'PRIORITY_D'
PRIORITY_4:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste
prioriteetidega
  shortValue:
    type: string
    description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:
      oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
  value:
    type: string
    description: Prioriteedi nimetus
  description: PriorityEnum 'PRIORITY_4'
PRIORITY_1A:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste
prioriteetidega
  shortValue:
    type: string
    description: Prioriteedi lühinimi
  department:
    type: array
    description: "Teenus(ed), kus kasutuses"
    items:

```

```
        oneOf:
        - type: object
        - type: boolean
        - type: string
        - type: number
    value:
        type: string
        description: Prioriteedi nimetus
    description: PriorityEnum 'PRIORITY_1A'
PRIORITY_INFO:
    type: object
    properties:
        orderNo:
            type: number
            description: Järjekorra number nimekirjades
        weight:
            type: number
            description: Prioriteedi kaal suhtes teiste
prioriteetidega
```

```
    shortValue:
        type: string
        description: Prioriteedi lühinimi
    department:
        type: array
        description: "Teenus(ed), kus kasutuses"
        items:
            oneOf:
            - type: object
            - type: boolean
            - type: string
            - type: number
    value:
        type: string
        description: Prioriteedi nimetus
    description: PriorityEnum 'PRIORITY_INFO'
PRIORITY_A:
    type: object
    properties:
        orderNo:
            type: number
            description: Järjekorra number nimekirjades
        weight:
            type: number
            description: Prioriteedi kaal suhtes teiste
prioriteetidega
```

```
    shortValue:
        type: string
        description: Prioriteedi lühinimi
    department:
        type: array
        description: "Teenus(ed), kus kasutuses"
        items:
            oneOf:
            - type: object
            - type: boolean
            - type: string
            - type: number
    value:
```

```

        type: string
        description: Prioriteedi nimetus
        description: PriorityEnum 'PRIORITY_A'
PRIORITY_1:
    type: object
    properties:
        orderNo:
            type: number
            description: Järjekorra number nimekirjades
        weight:
            type: number
            description: Prioriteedi kaal suhtes teiste
prioriteetidega
        shortValue:
            type: string
            description: Prioriteedi lühinimi
        department:
            type: array
            description: "Teenus(ed), kus kasutuses"
            items:
                oneOf:
                    - type: object
                    - type: boolean
                    - type: string
                    - type: number
        value:
            type: string
            description: Prioriteedi nimetus
            description: PriorityEnum 'PRIORITY_1'
PRIORITY_B:
    type: object
    properties:
        orderNo:
            type: number
            description: Järjekorra number nimekirjades
        weight:
            type: number
            description: Prioriteedi kaal suhtes teiste
prioriteetidega
        shortValue:
            type: string
            description: Prioriteedi lühinimi
        department:
            type: array
            description: "Teenus(ed), kus kasutuses"
            items:
                oneOf:
                    - type: object
                    - type: boolean
                    - type: string
                    - type: number
        value:
            type: string
            description: Prioriteedi nimetus
            description: PriorityEnum 'PRIORITY_B'
PRIORITY_2:
    type: object
    properties:

```

```

orderNo:
  type: number
  description: Järjekorra number nimekirjades
weight:
  type: number
  description: Prioriteedi kaal suhtes teiste
prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
  type: array
  description: "Teenus(ed), kus kasutuses"
  items:
    oneOf:
      - type: object
      - type: boolean
      - type: string
      - type: number
value:
  type: string
  description: Prioriteedi nimetus
  description: PriorityEnum 'PRIORITY_2'
PRIORITY_CONSULTATION:
  type: object
  properties:
    orderNo:
      type: number
      description: Järjekorra number nimekirjades
    weight:
      type: number
      description: Prioriteedi kaal suhtes teiste
prioriteetidega
shortValue:
  type: string
  description: Prioriteedi lühinimi
department:
  type: array
  description: "Teenus(ed), kus kasutuses"
  items:
    oneOf:
      - type: object
      - type: boolean
      - type: string
      - type: number
value:
  type: string
  description: Prioriteedi nimetus
  description: PriorityEnum 'PRIORITY_CONSULTATION'
  description: "PriorityEnums - PRIORITY_CONSULTATION,
PRIORITY_INFO, PRIORITY_1A, \
  \ PRIORITY_1, PRIORITY_2, PRIORITY_3, PRIORITY_4, PRIORITY_A,
PRIORITY_B, \
  \ PRIORITY_C, PRIORITY_D"
areaEnums:
  type: object
  properties:
    AREA_SOUTH:

```

```

        type: string
        description: AreaEnum 'AREA_SOUTH'
AREA_NORTH:
    type: string
    description: AreaEnum 'AREA_NORTH'
AREA_EAST:
    type: string
    description: AreaEnum 'AREA_EAST'
AREA_WEST:
    type: string
    description: AreaEnum 'AREA_WEST'
description: "AreaEnums - AREA_NORTH, AREA_SOUTH, AREA_EAST,
AREA_WEST"
deviceTypeEnums:
    type: object
    properties:
        DEVICE_TYPE_TETRA_PERSON:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_TETRA_PERSON'
        DEVICE_TYPE_MGIS:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_MGIS'
        DEVICE_TYPE_NUGIS:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_NUGIS'
        DEVICE_TYPE_MOBILE_CHIEF:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_MOBILE_CHIEF'
        DEVICE_TYPE_MOBILE:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_MOBILE'
        DEVICE_TYPE_EMAIL:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_EMAIL'
        DEVICE_TYPE_EMAIL_CHIEF:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_EMAIL_CHIEF'
        DEVICE_TYPE_DERIVED:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_DERIVED'
        DEVICE_TYPE_TETRA_VEHICLE:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_TETRA_VEHICLE'
        DEVICE_TYPE_TETRA:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_TETRA'
        DEVICE_TYPE_TETRA_CHIEF:
            type: string
            description: DeviceTypeEnum 'DEVICE_TYPE_TETRA_CHIEF'
description: "DeviceTypeEnums - DEVICE_TYPE_NUGIS,
DEVICE_TYPE_MGIS, DEVICE_TYPE_TETRA,\
    \ DEVICE_TYPE_TETRA_VEHICLE, DEVICE_TYPE_TETRA_CHIEF,
DEVICE_TYPE_TETRA_PERSON,\
    \ DEVICE_TYPE_MOBILE, DEVICE_TYPE_MOBILE_CHIEF,
DEVICE_TYPE_EMAIL, DEVICE_TYPE_EMAIL_CHIEF,\
    \ DEVICE_TYPE_DERIVED"

```